

Shop Order

Customer: City of Broken Arrow FD

Bryan Myrick

200 South 1st Street Broken Arrow, OK 74012 Dealer: Prof TX

Exp. Date: 06/07/2019 **Quote No:** 10101-0001

07/10/2019 13	18:24		Page 1
PART NO	S	DESCRIPTION	QTY
00-00-0010	ORDEF	R NUMBER	1
00-00-0106		Date: Post July 1, 2018	1
00-00-0113	Ambu	ulance Built to Star Of Life KKK-A-1882-F Certification, Latest	1
	Revision		
	Exceptions to	be noted on final sticker.	
00-10-0015		ntity (Vehicles)	1
00-91-0002	BODY NUM		1
00-91-0005	•	Y BODY DUE DATE :	1
00-91-0003		IN NUMBER:	1
00-91-0014		MANAGER: Jody Richardson	1
00-91-0027	_	E UNIT #: 170F22426XA	1
00-91-0028	RE-WRITE D		1
00-91-0040		POSITIONS MODS	1
00-91-0042		EAT BELT MOD:	1
00-91-0044		: 2 Positions	1
00-91-0050	•	SEAT BELT MOD:	1
00-91-0054	•	l: 4 Positions	1
00-99-9000		vel: 0 - Zero - ORIGINAL VERSION	1
00-99-9030	Change Orde	er Level: 0 - Zero - ORIGINAL VERSION	1
		Aluminum Composite Cabinet Configuration	1
01-1F-0000		ılar Ambulance	1
01-AO-1F48		F-550 4 x 2 DRW Reg Cab, 193" WB, XLT Turbo Diesel	1
01-1F-0085		: 6.7L Scorpion, Pwr Stroke V-8 Turbo Diesel - F-series	1
01-1F-0092		sis Front engine hood Insulation Not included by OEM Ford on all	1
	XLT		
	_	nge by Ford OEM Sept 2018,.	
	No Underhoo	od insulation installed by Ford OEM.	

07/10/2019 13:48:2		Page 2
PART NO S 01-1F-0101 01-1F-02CM 01-1F-02T5 01-1F-03PU 01-FL-TI06 01-FM-TI01 01-TU-0112 01-FM-EWH2 S	98R- Operator Command Regeneration Cab Interior Color: Medium Earth Gray F-series 2017 F-Series 193" wheelbase reg cab LWB GVWR: 18,000 pounds, FAWR: 7,000, RAWR: 13,660 (4x2) Tire, SPARE: Matching, Random Make Location: Shipped Loose Jack and Tire Tools: Install behind Passenger's seat Wheel Finish: Painted BLACK (6) Wheels ILOS PAINT: All (6) wheels shall be painted on BOTH sides, EXCLUDING the spare.	QTY 1 1 1 1 1 1 1 1 1 1
01-FM-DL01 01-FM-DTRL 01-FM-DTRQ 01-FM-TIR3 01-GD-0011 01-OE-47LF 01-OE-47A5 01-OE-47C1 01-RR-05B5 01-SY-0000 03-EA-47PR	Cab Door Locks: Power Daytime Running Lights: Ford OEM Ford F series High center mount brake light ordered with chassis Tires: F-4/550, All season, Random Mk 225/70 R19.5, LR F Ford Fleet 2019 Price level increase 47L - Ambulance Prep Package - F-Series Alternators Dual (OEM) Motorcraft Ford 67B 2017+ Ford F super duty Throttle: Intermotive High Idle Radio, F-series: OEM AM/FM/MP3 w/ 4 cab speakers Ford SYNC Basic installed Enhanced voice recognition Side Mirrors, OEM (F-Series), Pwr Glass, Heated, manual Telescopic XLT version mirrors Manually telescoping, folding trailer tow mirrors. Power adjustment of main section Heated glass, heated convex spotter mirror. Integrated clearance lamps with halo turn signals.	1 1 1 1 1 1 1
04-SU-1101 04-SU-1729	Ford OEM Front Suspension Coil Spring Std F-450/550 Rear Suspension: F-550 Liquid Spring, 2017+ *NOTE: If Oxygen Lift is installed, then if the O2 Cmpt Door is open, the kneeling suspension will not dump until door is closed* Part; SUSP-DS147FS3 Rated for 13,660 GAWR	1
01-TT-0208	""DETAIL"" Liquid Rear Suspension Decal Install Prior to delivery, detail department is to install a Black Decal on the dash:	1
04-SU-0601 04-SU-0610 04-SU-0651 04-SU-03F6 ZZ-ZZ-ZZZZ	Kneeling Feature: Enable Switch Located in CAB console Kneeling Feature: Activated by TRAILING rear access door Exhaust system termination point: OEM Location, Rt Rear Rear Stabilize (Anti-Sway) Bar: OEM, Ford F-450/550 == 172 x 96 T-1 LWB Ford X-Series, Conversion - 15.009 05/07/19 == BODY Generation 9	1 1 1 1 1
01-X0-0001	Body Certification Provided to Specification The Body Manufacturing group certifies that this Ambulance Modular body has been constructed of materials and methods that meet or exceed the standards as required by the AEV Generation 9 Body Specification of December 11th, 2018 as published. The body shall be built and constructed to the submitted work order and the AEV Generation 9 build standards without exception. Any Deviations to the work order or Generation 9 build standards must have written authorization	1 1

07/10/2019 13			Page 3
PART NO	S	prior to body delivery to AEV by the AEV Body Engineer on file.	QTY
02-00-0002	XS		1
02-00-0003		ALL DIMENSIONS ON ORDER WILL HAVE A VARIABLE OF + OR - 1/4" Unless specified as a minimum dimension.	1
02-00-0010 02-B2-X092 02-00-0006 02-00-0060 02-00-0010 02-B0-CC01	S	Body, Mod: Type 1 LWB, 174" x 96" x 72" Interior HR - 6" Body Drop x Headroom: 74/72" Finished Interior - Ducted A/C Compartment and Entry door Full length Stainless Steel Hinges std Compartment Construction: STANDARD, Unless Specified Otherwise Sides and Back: .100 Polished Aluminum Diamond Plate Compartment Ceiling: .090 Aluminum Sheet	1 1 1 1 1
02-B0-SFL0		Compt Floors: Sweep-out, even with bottom door jamb Unless Specified Otherwise.	1
02-B0-CC41 02-B0-CC18 02-B0-CEX2 02-B0-09B2 04-EA-09B1 02-B0-CESL	S	Compartment floor thickness .125" standard Compartment Ventilation - Austin Flap style #VT-2495-A01, STD Curbside Entry Door (CSE): 81.813"H x 31.0"W X w/6" drop C/S Door Check : Double Action Gas Shock Door Swing Angle: Set to 90-100 degrees Curbside Entry Door (CSE) Location: Standard, Rearward of M-7 Compt	1 1 1 1 1
ZZ-ZZ-ZZZZ 02-B1-SWX4 05-IL-09X1		if Optioned Step Well, CURBSIDE Entry, 2-Step Embossed D/P, Elongated for Portable O2 w/door Light, Step Well: 3" Kinequip #4302S, Chrome, LED, STD Light shall come on with CS door only.	1 1
02-BC-0604 02-BC-0700		Talk through, Cab to Mod Window (T1) With Sliding Window 2017+ F Cab Roof Support : 3/16" x 3" Steel Plate, Per Engineering Drawing Between cab headliner and bottom side of cab roof to prevent oil canning.	1
02-BD-0120 02-X1-M149	S	Body Drop: 6" Both Sides, Ahead of Rear Wheels M-1 Compt (LF): 88.8" H x 25" W x 21.0" D NOTE: Backboards on the right hand side. There shall be a strap provided with buckle to prevent backboards from falling out when door is opened. The left side of this compartment shall have at the floor height Transverse O2. Above cylinder storage area storage for a Stryker Stair-Pro model# 6252. A strap with buckle to prevent stair chair from falling when door is opened.	1 1
02-B0-CCM1		/ M-1 Compartment Construction: Smooth Aluminum ILOS Sides and Back: .090 Aluminum Sheet Compartment Ceiling: .090 Aluminum Sheet	1

07/10/2019 13:		25	Page 4
PART NO	S	DESCRIPTION	QTY
02-B0-CC09		Compartment Door Panel: Smooth Aluminum	1
02-B0-CFM1		/ Compartment Surface Finish: Speed Liner, Black - M-1	1
02-B0-CFN1		Component Finish (Shelf/Tray/Divider): Speedliner Black Each	2
02-B0-CC16		Compartment Door Ventilation - small punched half moon Louvered Door	1
02-B0-SWFL		Compt Floor: Sweep-out, even with bottom door jamb.	1
02-B0-CC41		Compartment floor thickness .125" standard	1
02-B0-CC41 02-BC-10A5		· · · · · · · · · · · · · · · · · · ·	' '
		Compartment Door: SINGLE DOOR, Forward hinged 2-point Latch	1
02-B0-09B0		Compt Door Check: Double Action Gas Shock	1
04-EA-09B1		Door Swing Angle: Set to 90-100 degrees	1
04-TS-0040		Shelf requested In M-1 oxygen Compartment - Exception to 3.11.3	1
		required	
04-TS-0042		Divider requested In M-1 oxygen Compartment - Exception to	1
		3.11.3 required	
04-TS-1400	S	Divider, Vertical, Full compartment height, alum, Fixed	1
		NOTE: Installed next to Transverse O2 for Backboard Storage to Rear of	
		Divider.	
		Make Sure the Divider is FULL DEPTH	
04-TS-STR1	S	Strap: 2" Webb, w/ Chrome Metal Seatbelt buckle footman	1
04 10 01111	0	loops	′
		This strap is to secure things from moving while the vehicle is in motion, it	
		is not a tested device and is not designed to be a crash stable securing	
		device.	
		NOTE I A II I E WALLET BY II A GOLD OF	
		NOTE: Installed From Wall #1 to Divider for Stair Chair Storage	
04-TS-STR2		2" Footman Loops: Locate back from jamb	1
11-X0-0026		Stowage label - 2 inch Strap Seatbelt Matl /Seatbelt buckle	1
		w/loops 20lbs	
04-TS-STR3	S	Strap: 2" Webb, w/ Chrome Metal Seatbelt buckle w footman	1
		loops IATS	
		This strap is to secure things from moving while the vehicle is in motion, it	
		is not a tested device and is not designed to be a crash stable securing	
		device.	
		NOTE: From Divider to Wall #3 for Backboard Storage	
		NOTE. I folli bivider to wall #3 for backboard Storage	
OA TO OTDO		1 21 Factman Lagran Lagran back from James	
04-TS-STR2		2" Footman Loops: Locate back from jamb	1
11-X0-0026		Stowage label - 2 inch Strap Seatbelt Matl /Seatbelt buckle	1
		w/loops 20lbs	
11-YZ-0903		Non-Compliant to 3.11.3 KKK-F Storage in Oxygen	1
		Compartments Notice	
04-TS-11D0	S	Shelf, Fixed, Ext: .125 Aluminum, 2" Upward Lips	1
		NOTE: Install Fixed Shelf From Wall #1 to Divider Above Transverse O2 for	
		Stair Chair Storage.	
		Shelf to be 17" W x 21" D (Previous Shelf wasn't deep enough for stair chair	
		to sit on)	
		,	
11-YZ-0903		Non-Compliant to 3.11.3 KKK-F Storage in Oxygen Compartments	1
1. 12 3333		Notice	′
i			

07/10/2019 13:			Page 5
PART NO 05-EL-40X1	S	Light, Compt, LED, M-1 (LF): Intertek Model No ZY-156-LEDC, STD MTB Note: 1" Hole, Centered in Compartment ceiling	QTY
08-B1-02LI 08-B1-02LI	X	Prep for 02 Lift: NONE Transverse Oxygen Compartment in M1 Near Wall #1 Above Floor Line	1
09-ZZ-ZZZZ 02-X1-M275 02-B0-CCM2	S	M-2 Compt (LFM): 34.5" H x 48" W x 20.0" D M-2 Compartment Construction: Smooth Aluminum ILOS Sides and Back: .090 Aluminum Sheet Compartment Bottom: .125 Aluminum Sheet Compartment Ceiling: .090 Aluminum Sheet	1 1 1
02-B0-CC09 02-B0-CFM2 02-B0-CFN1 02-B0-SWFL 02-B0-CC41 02-B1-BION ZZ-ZZ-ZZZZ		Compartment Door Panel: Smooth Aluminum Compartment Surface Finish: Speed Liner, Black - M-2 Component Finish (Shelf/Tray/Divider): Speedliner Black Each Compt Floor: Sweep-out, even with bottom door jamb. Compartment floor thickness .125" standard Bio-Waste None;	1 1 2 1 1 1
	S	Electrical Compartment : Circuit Board Storage relocated to M 2.5 Compartment Door hinged at bottom.	1
02-BC-10D0 02-B0-09B0 04-EA-09B1 04-TS-11D3	S	Doors, Compartment, DOUBLE DOORS (std) Compt Door Check: Double Action Gas Shock Door Swing Angle: Set to 90-100 degrees M-2 Shelf, Adjustable, Ext: .125 Alum, 2" Upward lip	1 1 1 1
04-TS-4010 04-TS-2100 04-TS-21A2	S	Shelf Bracket: CPI 90 Deg with Slotted Holes Tray, Roll-out: Std Duty .125 Thick tray, Fixed position Slides: Pair of Accuride AH760 250 lb 18" Full Extension with Right-side locking One locking slide with yellow thumb catch release (1) SLID-9301Non-locking (1) SLID-9308 Right side locking Right sias typical 18"	1 1
05-EL-40X1		/ Light, Compt, LED, M-1 (LF): Intertek Model No ZY-156-LEDC, STD MTB Note: 1" Hole, Centered in Compartment ceiling	1
05-HA-13HA 05-HA-13I0 04-HA-1310	S	/ / Condenser, 12V: ACC #X717, Dual Fan, Front of Module / ARCTIC WEDGE Condenser Cover: Install on Front Body over Condenser part #WEDGE-HVAC Use of the wedge is subject to engineering review and approval for Light duty Type III models. A reduction in light sizing above the wedge may be required.	1 1 1
12-B0-0004 06-RR-13Z0 06-RR-13O4	s	Painted Part- Painted Special color to match Body Standard Conduit: 1-1/2"", with pull wire CONDUIT ORIGINATION POINT: M-2.5 compartment	1 1 1

07/10/2019 13:			Page 6
PART NO	S	DESCRIPTION	QTY
06-RR-13T3		CONDUIT TERMINATION POINT: Behind A/A Board (Panel)	1
06-RR-13Z1		Conduit #2 1-1/2" diameter, with pull wire, Type 1 units Std	1
06-RR-13O1		CONDUIT ORIGINATION POINT: Electrical Circuit board cabinet	1
06-RR-13T2		CONDUIT TERMINATION POINT: Behind Driver's seat	1
09-ZZ-ZZZZ		[*]	1
02-X1-M295		M-2.5 Exterior electrical (LWW): 28" H x 48" W x 6" D	1
02-B0-CC03		/ Compartment Construction: 100% Smooth Aluminum ILOS	1
02 20 0000		Sides and Back: .090 Aluminum Sheet	-
		Compartment Bottom: .125 Aluminum Sheet	
		Compartment Ceiling: .090 Aluminum Sheet	
		Compartment Cening090 Adminum Sheet	
12 DC CA20	C	L Compositment Finish: DA Cond	4
12-DC-GA20	3	Compartment Finish: DA Sand	1
		Paint Shop to Refer to Body Section of Shop Order for Details.	
02 BC 40C4		L Door Comportment SINCLE DOOD 2 noint Lotab II OS Ton Hingard	4
02-BC-10C1		Door, Compartment, SINGLE DOOR, 2-point Latch ILOS Top Hinged	7
		Gas strut hold opens on both sides to lift door.	
		Nylon strap installed on foot man loop to aid in closing door after opened.	
02-B0-09B0		Compt Door Check: Double Action Gas Shock	1
04-EA-09B1		/ Door Swing Angle: Set to 90-100 degrees	1
05-EL-40X2	S		1
00 1011_		MTB Note: No light hole required	
		(2) on ceiling, near door jamb, on angled bracket	
		(2) on coming, near accir james, on angiou bracket	
09-ZZ-ZZZZ		I	1
02-X1-M375	S	1	1
02-X 1-W373	3	NOTE: Duo Therm TO BE ABOVE M-3	′
		NOTE. Duo Tiletti TO DE ADOVE M-3	
		NOTE: Box Off in M3 Compartment for Interior Glove Storage	
		(Make Close off around Gloves As Small as Possible - SEE PHOTO ON FILE)	
		(Make close of around cloves As official as 1 ossible off 11010 off 1122)	
		To house SCBA Bracket and Vulcan Flashlight Below Shelf	
		To house GodA Bracket and Valean Flashinght Below Offich	
02 00 00 00		M 2 Compartment Construction, Smooth Aliminum II OS	4
02-B0-CCM3		M-3 Compartment Construction: Smooth Aluminum ILOS	1
		Sides and Back: .090 Aluminum Sheet	
		Compartment Bottom: .125 Aluminum Sheet	
		Compartment Ceiling: .090 Aluminum Sheet	
02-B0-CC09		Compartment Door Panel: Smooth Aluminum	1
02-B0-CFM3		/ Compartment Surface Finish: Speed Liner, Black - M-3	1
02-B0-CFN1		Component Finish (Shelf/Tray/Divider): Speedliner Black Each	1
02-B0-SWFL		Compt Floor: Sweep-out, even with bottom door jamb.	1
02-B0-CC41		Compartment floor thickness .125" standard	1
02-BC-1000		Compartment Doors: DOUBLE DOORS, ILOS single	1
			1
02-B0-09B0		Compt Door Check: Double Action Gas Shock	7
04-EA-09B1	_	Door Swing Angle: Set to 90-100 degrees	1
04-TS-22M3	S	Box Off For Interior Glove Cabinet	1
04-TS-22M3	X	Corner Bracket for SCBA Bracket to Mount	1
		Install Corner Bracket on Wall #1 and Wall #2 for SCBA Bracket to Mount to.	

DESCRIPTION 05-EL-40X3 Light, Compt, LED, M-3 (LR): Intertek Model No ZY-156-LEDC, STD MTB Note: 1" Hole, Centered in Compartment ceiling M-5 Compt (RR): 28" H x 30" W x 20" D NOTE: to House SBA Bracket and Vulcan Flashlight M-5 Compartment Construction: Smooth Aluminum ILOS	1 1 1
09-ZZ-ZZZZ 02-X1-M574 S M-5 Compt (RR): 28" H x 30" W x 20" D NOTE: to House SBA Bracket and Vulcan Flashlight	1 1
02-X1-M574 S / M-5 Compt (RR): 28" H x 30" W x 20" D NOTE: to House SBA Bracket and Vulcan Flashlight	1 1
NOTE: to House SBA Bracket and Vulcan Flashlight	1
	1
TUZ-DU-GGNO I NI-O GOMOANMENI GONSINIGNON: SMOOM AUMINUM ILUS	
Sides and Back: .125 Aluminum Sheet Compartment Bottom: .125 Aluminum Sheet	
Compartment Ceiling: .090 Aluminum Sheet	
02-B0-CC09 Compartment Door Panel: Smooth Aluminum 02-B0-CFM5 Compartment Surface Finish: Speed Liner, Black - M-5	1
02-B0-CFN0 Component Finish (Shelf/Tray/Divider): Standard	1
02-B0-CC12	1
02-B0-CC16 Compartment Door Ventilation - small punched half moon Louvered Door	1
02-B0-SWFL Compt Floor: Sweep-out, even with bottom door jamb.	1
02-B0-SWFL X Corner Bracket for SCBA Bracket to Mount	1
Install Corner Bracket on Wall #2 and Wall #3 for SCBA Bracket to Mount to.	
02-BC-1000 Compartment Doors: DOUBLE DOORS, ILOS single 02-B0-09B0 Compt Door Check: Double Action Gas Shock	1
05-EL-40X5 Light, Compt, LED, M-5 (RR): Intertek Model No ZY-156-LEDC, STD MTB Note: 1" Hole, Centered in Compartment ceiling	1
09-ZZ-ZZZZ	1 1
02-X1-M775	1
02-B0-09B0 Compt Door Check: Double Action Gas Shock 04-EA-09B1 Door Swing Angle: Set to 90-100 degrees	1
05-EL-40M7	1
02-X1-M846 S / M-8 Compt (RF): Batt Compartment 17" H x 25.25" W x 20" D	1
NOTE: To house Inverter, Suction Pump, Breaker Box, Etc.	
02-B0-CCM8	1
Compartment Ceiling: .090 Aluminum Sheet	
02-B0-CFM8 Compartment Surface Finish: Speed Liner, Black - M-8 02-B0-CC16 Compartment Door Ventilation - small punched half moon	1
Louvered Door 02-BC-10A0 Compartment Door: SINGLE DOOR, 1-point Latch	1

07/10/2019 13:			Page 8
PART NO	S	DESCRIPTION	QTY
02-B0-09B0		Compt Door Check: Double Action Gas Shock	1
04-EA-09B1	_	Door Swing Angle: Set to 90-100 degrees	1
05-EL-40X5	S	/ Light, Compt, LED, M-8 (RF): Intertek Model No ZY-156-LEDC, STD MTB Note: No Light Hole Required	7
		NOTE: Mount Light on Wall #3	
<i>ZZ-ZZ-ZZZZ</i>		!	1
02-X1-RA02		Rear Access Doors: 46.75" Wide x 60.625" High	1
09-IN-3STX		Insulation: All Module Doors, 2" Black Foam, X-Series units #INSU-ARM2	1
09-IN-3STY		Insulation: All Module Doors, Black Sound Proofing, X-Series units #INSU-KSU	1
ZZ-ZZ-ZZZZ			1
		WELDING - CELL 1	1
04-AS-06A0	_	Fuel Fill: CPI C1045, Open Housing, Polished Bezel - (Std)	1
04-AS-0200	S		1
		12" Wide diamond plate shield shall be contoured to fit under the fuel filler housing - ALL THE WAY TO THE SKIRT RAILS	
		Housing - ALE THE WAT TO THE ORIGINALES	
04-AS-0801		Urea Fill: CPI FG5301-1 Inlet Housing	1
04-AS-0420	S		1
		9-1/2" Wide diamond plate shield shall be contoured to fit under the filler	
		housing ALL THE WAY TO THE SKIRT RAILS	
04-AS-0370		Wire/Hose Cover : Diamond Plate, Between Cab & Module	1
ZZ-ZZ-ZZZZ			1
04-AS-0382		Pass-Thru cab mounted wind deflector - Aluminum Diamond Plate	1
		AEV Weld shop builds and installs, Attaches to Cab side of bellows mount. Same width as bellows mount with short down turned sides. Secured on sides to	
		bracket. Small up-turned lip at back of wind deflector.	
		non-Painted	
04-BW-ABY1		Bumper, Rear: X-Series 96" Wide w/ Embossed D/P Pontoon Covers	1
		X-Series Design. Purchasing NOTE: Special Order Delete standard rear bumper assembly from	
		Mickey Body	
04-BW-FLIW		Step, Center:2" x 9" Bar Grate, flip-up, recessed	1
		Standard with X series, CCP purchased item	
		MTB delete flip up step.	
04 EA 4206		L Dook Bumperer Bubber D Chang All Lly 40" Wy 4" Thick V against	
04-EA-1306 04-BW-RF20		Dock Bumpers: Rubber D-Shape, 4" H x 18" W x 4" Thick X -series Fenders, Rear: Rubber Extruded, Black	1
57 DW-1(1 20		MTB #7215439	'
04-BW-DN00		Corner Caps "X Series:	1
04-BW-DN23		Front Corner Caps: 24.0" High, Alum Diamond Plate	1
		NOTE: Standard Front Corner Caps from Mickey Body	

07/10/2019 13:48:2		Page 9
PART NO S	DESCRIPTION	QTY
04-BW-DN33	Rear Corner Caps: 24" High, Polished Diamond Plate, X-Series Purchasing NOTE: Special Order from CCP. Delete standard rear corner caps from MTB	1
ZZ-ZZ-ZZZZ 04-BW-DP2A XS	Rear Kick Plate: Smooth Aluminum, Recessed, X-Series ILOS for Chevrons Purchasing NOTE: Special Order from CCP. Delete standard rear kick-plate from Mickey Body	1 1
04-EA-2500	/ Rear Recovery Eyes: (2) Cast Iron 1 3/16" center eye Paint Black 1-3/16" ID center eye	1
04-EA-2501 04-BW-TA03 S 04-BW-TA04 05-EL-44TO ZZ-ZZ-ZZZZ	Eye Access Holes: (2) 5 inch Dia, Align w/ eyes Recessed Tag Area: Smoot Aluminum Location: Centered in the kick plate Tag Lights: (2) Kinequip LED #132703C, ILOS	1 1 1 1
05-EL-44XL	Rear Flip up Step Lights, (2), Kinequip LED Strip light #KFL-SLC8 Mounted above rear bumper flip up step, in kick-plate recess	1
ZZ-ZZ-ZZZZ 04-BW-DP93	 Front Stone Guards: 24.0" High, Alum Diamond Plate If there is a body drop, add the appropriate height to Stone Guards.	1
04-BW-DPX5	/ Skirt Rails: Polished D/P, W/ Whelen Strip Lights, "X Series" ILOS X-Series Design	1
	Purchasing NOTE : Special Order from CCP. Delete standard skirt rails from Mickey Body	
05-EL-SKWT	<i>J</i> (6) Red Skirt rail LED Strip Lights, Whelen PSR00FRR, Red Lens Lite-PSR00FRR strip light Red LED, Red LENS. Any lights placed forward of Rear axle to not steady burn due to NHTSA and	1
	FMVSS requirements.	
05-FS-10P4 11-ZZ-0909	/ Flash Pattern: None, Internal Light Head Flasher / CAAS GVS Exception to C.8.2.12 Warning lights non-compliant to Table 1 dual mode The Warning lights pattern chosen by the Customer does not meet either the required photometric output flash or the pattern required of Table 1, Dual mode, Primary Secondary of CAAS C.8.2.12. An exception has been taken for only this	1
05-EL-SKX7	 Skirt rail Lights, Front Pair (1) Each side, forward position, in front skirt rail 	1
05-EL-SKX8	/ Skirt rail Lights, Middle Pair (1) Each side, rear position, in front skirt rail	1

07/10/2019 13:			Page 10
PART NO 05-EL-SKX9	S	J Skirt rail Lights, Rear Pair (1) Each side, in rear skirt rail	QTY 1
05-EL-SKZ1 05-EL-SLGM 05-EL-SLGZ	S	Skirt rail Lights Switched: Separate Switch in the Cab console LED Ground Lights: (6) Whelen White 8" LED lights, In Skirt rails LED Ground Lights Switched: Separate Switch in the Cab console and with any entry or compartment door.	1 1 1
05-IL-09S0	S	/ Ground lights- Whelen PSC0CDCR LITE-PSC0CDCR strip led light LED	6
04-DE-0010 04-EA-09A0 04-EA-14X5		Rear Deflector: NONE (2) Rear Door Hold Opens: Grabber Style, each door Running Boards:Embossed D/P,Bar-Grate Insert, X Series - 2 Door X- STD F 2017 Bar grate insert 16" long BOAR-XF2017	1 1
04-EA-1122 04-EA-1132 05-EL-39X4		Mud Flaps Front: Modular, Rubber with AEV/Rev Logo std Mud Flaps Rear: Modular, Rubber AEV LOGO Lights, Cab Step: 2" White LED, surface mount on box front, (#TOCACCCR) - x ser Chrome Flange Activate by Cab Doors. Constant Hot. (2) LITE-TOC (2) Lite-TFLANGEC	1 1 1
ZZ-ZZ-ZZZZ 05-EA-2110 05-EL-1450 05-EL-1201 05-PH-LF20 05-EL-1800 05-EL-18F0 05-EL-18FX	S	EXTERIOR LIGHTING / HEAT/AC / INSULATION - CELL 2 Grill Guard: Ranch Hand - Model TBD Front Turn, Whelen M6 Series w/arrow, LED, Pair, IATS Location: Front of Module, Below Light Bar M6 Series Flange, Chrome Plated, Std: Each CORNER CAP LED ICC/WARNING LIGHTS Front Corner Cap LED ICC/Warning Lights: Warnings RED/WHITE Front Center ICC Lts: (3) AMBER Kinequip LED #KML-A-SC Switch w/ Headlights.	1 1 1 1 2 1 1
ZZ-ZZ-ZZZZ 05-EL-18R5 05-EL-18RX		 Rear Corner Cap LED ICC/Warning Lights: Warnings RED/AMBER Rear Center ICC Lts: (3) RED Kinequip LED #KML-R-SC Switch w/ Headlights	1 1 1
ZZ-ZZ-ZZZZ 06-SW-CC01 ZZ-ZZ-ZZZZ 05-EL-2503	S	ICC Warning Lights Switched: Primary Only LED Lights, Patient on Board, X-Series Units Install lights in the following locations: M-1 Compartment - On Door panel, Outer Edge, Centered Top to Bottom	1 1 1 1
		M-2 Compartment - On Forward Door panel, Upper outer corner M-2 Compartment - On Trailing Door exterior, Upper outer corner	

07/10/2019 13:48	5.27	Page 11
PART NO S		QTY
	M-3 Compartment - On Forward Door panel, Outer Edge, Centered Top to	
	Bottom	
	M-3 Compartment - On Trailing Door exterior, Outer Edge Centered Top to	
	Bottom M.F. Compostment On Forward Deer name! Outer Edge Univer Outer Corner	
	M-5 Compartment - On Forward Door panel, Outer Edge, Upper Outer Corner M-5 Compartment - OnTrailing Door, Outer Edge Upper Outer Corner	
	M-7 Compartment - On Door panel, Outer Edge, Centered Top to Bottom	
	M-8 Compartment - On Door Panel Outer Edge Upper Outer Corner	
	C/S Entry Door - Below Center panel, Outer Edge	
	Rear Entry Doors - (1) each door, Below Center panel, Outer Edge	
	A Total of (12) Lights	
	A Total of (12) Lights	
05-EL-252C S	S Red/White LED POB Light, Kinequip Each, X-Series Units	12
00 LL 2020	NOTE: POB Installed on Compartment Doors. M1 Door Panel Centered on Outer Edge, M2 Front Door Panel Upper Outer Corner, M2 Rear Door Exterior Upper Outer Corner, M3 Front Door Panel Above Door Handle Height, M3 Rear Door Exterior Above Door Handle Height, M5 Front Door Panel Upper Outer Corner, M5 Door Exterior Upper Outer Corner, M7 Door Panel Below Handle Height and M8 Door Panel upper Outer Corner. One Installed on Each Entry Door Panel Below Center Handle Panel Toward Outer Edge. Total	72
05-EL-256X 05-EL-SKWR S	of (12) Lights Switched with Entry Door Switch (2) Red Skirt rail LED Strip Lights, Whelen PSR00FRR, Red Lens Lite-PSR00FRR strip light Red LED, Red LENS.	1 1
	NOTE: Install on Cab Doors to Flash When Doors are Open	
05-FS-10P4 11-ZZ-0909	Flash Pattern: None, Internal Light Head Flasher CAAS GVS Exception to C.8.2.12 Warning lights non-compliant to	1
	Table 1 dual mode	
	The Warning lights pattern chosen by the Customer does not meet either the required photometric output flash or the pattern required of Table 1, Dual mode, Primary Secondary of CAAS C.8.2.12. An exception has been taken for only this element to the standard.	
ZZ-ZZ-ZZZZ	I	1
05-EL-42XX	TAIL LIGHTS	
05-EL-23X2 05-EL-35A6	Tail Lights, Whelen LED Combinations , STD, X-Series Side Marker Lts: (2) Whelen M6 LED Red Running/Turn, ILOS	1 1
	Running to be wired as Tail and Turn to be wired as Stop. These are to both have flanges.	
05-PH-LF20 05-EL-2414 05-EL-2420	M6 Series Flange, Chrome Plated, Std: Each Whelen M6 Series, Tail Lights Stop/Tail, Whelen M6 Series, LED, Pair	2 1 1

07/10/2019 13:	.40.4	21	Page 12
PART NO	S	DESCRIPTION	QTY
05-EL-2426		Turn, Whelen M6 Series, LED, Pair	1
05-EL-2432		/ Back up, Whelen M6 Series, LED, Pair	1
05-PH-LF24		M6 Series Flange, Chrome Plated, Std: Three Pair	1
06-EC-4326		Note: Standard Back-up Alarm Installed	1
		•	1
06-EC-4331		Access Panel: NONE Required	7
ZZ-ZZ-ZZZZ			1
05-EL-4300		Patient Area Turn and Stop Light Indicator: None	1
05-EL-43SU		Third (3rd) Brake Light: Whelen RED LED Strip Light #PSROOXRR	1
05-EL-4584	S	SCENE/FLOOD LIGHTS	1
05-EL-46DL		Rear Load Lights: Built into Rear Light Bar	1
06-EL-18RE			1
	0	Activate: Rear Flood Switch, Reverse and Lead RA Door	1
06-EL-1800		Add Activation: All Scene and Load lights. come on w/ Reverse	7
06-EL-SST5	S		1
		chrome flange	
		Switches in console labeled left and right flood.	
		(Discontinued June 2018)	
		(2.000.1111.000.00110.2010)	
05 DU L T55		La cotto m (4) Contour de mande de la face de la	4
05-PH-LT5F		Location:(1) Centered on each side of module	7
06-EL-SSS7		Switches Labeled: (1) Left Flood & (1) Right Flood	1
ZZ-ZZ-ZZZZ		 	1
05-HA-0650		Secondary Heat/AC: Duo-Therm 41001-611 CoolCat Heat Pump	1
		Install condenser air shroud.	
05-HA-0660	S	Secondary AC Location: Over M3 compartment	1
05-HA-13XG		HVAC SYSTEM: ACC Hot Weather Pkg, 2nd Comp, Ducted - Above ALS -	1
03-11A-13AG	3		,
		X Ford ILOS	
		Unit is at the I-1 position of the ALS	
		Includes 2nd AC compressor installed by ACC	
		NO Heppa MAX system	
		electric water valves included.	
05-HA-0406		Water Pump: Installed required	1
05-HA-1306		A/C Condenser: See M2 Cmpt for Condenser Information	1
			1
05-HA-14A2		Return Air Plenum: Walk Thru Side of RF-ALS	1
		Stamp Powder Coated Steel Vent (VENT-FLOORV) mounted near floor. NOTE:	
		Recessed pocket door or Walk-thru side cabinets not available.	
05-HA-14A0		Side Plenum Grille, Return Air: Stamped Powder Coated Steel	1
05-HA-14C8		/ AC Compressor, 2nd AC installed 2017 Ford F-450/550 included	1
		Vehicle transported to facility for 2nd AC compressor, Bracket, Drive belt and	-
		installation. Coordinate with engineering at earliest opportunity prior to production	
		assembly start	
05-HA-14P1		AC Evaporator: ACC (AC/Heat unit) w/ dual fans - Hot Weather	1
		Horizontal	
05-HA-1404		Heater Hoses: EPDM - Nomex Rubber (per Ford QVM)	1
05-HA-14A0		/ Side Plenum Grille, Return Air: Stamped Powder Coated Steel	1
05-HA-14B1		Filter, Washable Carbon Pre-Filter	1
05-HA-15A0		Ther, Washable Carbon Tre-Ther Ducted AC Delivery: insulated & foil wrapped, 10 registers	1
	v		1
05-HA-1801	X	Contral HVAC Thru LX1 - See Quincy for Details	7
05-HA-1820		Air Curtain: Side Door - NONE	1
10-HA-0304		Thermostat, Rear A/C, Digital, Class 1 LX-1, Includes Front Also	1
		Thermostat located in A/A as part of LX1 Multiplexed System	

07/10/2019 13	:48:	27	Page 13
PART NO	S	DESCRIPTION	QTY
05-IL-00X1		Ceiling Panels; white Aluminum Composite (2) piece X series MATL-CEILING	1
05-IL-01XR		Dome Lts, LED "LIGHT RAIL" system MTB NOTE: Requires X-Series Roof Structure design	1
05-IL-0396 05-IL-05L1 ZZ-ZZ-ZZZZ	s	(2) Lights, 12VDC, Thinlite #716XL, 18" Fluorescent, Recessed LOCATION: Curbside Ceiling Panel 	1 1
05-IL-05T4	S	/ Check out Light Switch: Momentary For Multi-plex electrical systems.	1
		Part SWIT-C05 (1) Momentary switch to be installed on Curbside Wall	
06-EC-0916		<i>J 30 Minute Activation</i> Light Circuit to stay on for 30 minutes when entry door is opened.	1
06-EC-09P1		Timer to Power: Street side DOMES, High intensity	1
77-77-777 06-EC-2500	C		1
06-EC-2500 06-EC-09P1	S	Timer: Through LX1 System (No Switch) Timer to Power: Light Rail, High intensity	1
ZZ-ZZ-ZZZ			1
06-RR-0701		(2) Speakers: Stereo, patient compartment, in pad over rear doors	1
10-RR-0600 ZZ-ZZ-ZZZZ		Volume Control, Stereo: in A/A, For Rear speakers 	1
09-MH-0801		IV Hook No 1: CPI Rubber IV 2008-1(Recessed Mount), ILOS	1
29-MH-1000		/ LOCATION: Over head/chest area, primary patient on COT	1
09-MH-0802		IV Hook No 2: CPI Rubber IV 2008-1(Recessed Mount), ILOS	1
09-MH-08P2 ZZ-ZZ-ZZZZ		LOCATION: Over head/chest area, secondary patient on S/B 	1
09-MH-0803		IV Hook No 3: CPI Rubber IV 2008-1(Recessed Mount), IATS	1
09-MH-08P3		LOCATION: Over Knee/Waist area, primary patient on COT	1
09-MH-0804 09-MH-08P4		IV Hook No 4: CPI Rubber IV 2008-1(Recessed Mount), IATS LOCATION: Over Knee/Waist area, secondary patient on S/B	1
09-MH-2AC6		Recessed C/S Grab Rail, ceiling: 1.25 Dia 3 pt, 72in, Yellow	1
		Antimicrobial	
		Grab Rail will be recessed in a ABS pan.	
09-MH-2AYX		/ Recessed S/S Grab Rail, ceiling: 1.25 Dia 3 pt, 72in, Yellow	1
		Antimicrobial	
		Grab Rail will be recessed in a ABS pan.	
		Insulation DVO: V Codes with	
05-IN-3STX		Insulation PKG: X-Series units Insulation: All Module Doors, 2" Black Foam, X-Series units #INSU-ARM2	1 1
05-IN-3STY		Insulation: All Module Doors, Black Sound Proofing, X-Series units #INSU-KSU	1

07/10/2019 13:		21	Page 14
PART NO 05-IN-4ST9	S	Insulation Sound Deadening: Generation 9 Floor Generation 9 Floor Design	QTY 1
05-IN-5STX		Insulation: Walls/Ceiling, 1" Black Foam, X-Series units #INSU-ARM1	1
05-IN-6STX		Stepwell Insulation: Urethane Froth Insulation, STD X-Series units MTB to spray urethane froth insulation, underside of step well area	1
		MTB option #7208738	
05-IN-6SXS		Black Sound Proofing, Doors & Select Walls, X-Series units INSU-KSU	1
05-IN-6SXT		Insulation: Circumferential PKG, 1LEF Reflective, X-Series units #KSU-F1250	1
05-PA-LB0F 05-PA-LB64	S S	Front Light Bar: Whelen, 4500 Light Bar With Rota Beams, 86" L Light Bar LOCATION: Box FRONT / Surface mount Body skin to be reinforced by 1/4" plate secured to tubes in at minimum two locations.	1 1
		NOTE: Center Rotabeam fixtures set to fast rotate (175 rpm) the others will be standard rotate (95 rpm)	
05-PA-LB6Y	S	/ Front lightbar flash pattern please define Center Rotabeam fixtures set to fast rotate (#8 Rotator 150 CCW) the others will be standard rotate (#4 Rotator 75 CCW)	1
		The front of the modular body shall be equipped with a Whelen 45B88 Rota Beam light bar to include Four [4] 45ADRBR, Four [4] 45ADRBB and Two [2] 45ADRBC.	
06-SW-0006 06-SW-0007 06-SW-0008 06-SW-0011 06-SW-0011 05-PA-LB0K 05-PA-LB62	X S S	/ Light Bar LOCATION: Box REAR / Surface mount Body skin to be reinforced by 1/4" plate secured to tubes in at minimum two locations.	1 1 1 1 1
		LIGHT NOTE: Whelen 45B88 Rota Beam light bar to include Four [4] 45ADRBR, Two [2] 45ADRBB, Two [2] 45ADL8A and one [2] 45LRLOAD	
05-PA-LB6Z 06-SW-0006 06-SW-0007 06-SW-0008 06-SW-0011		Rear lightbar flash pattern please define Switches Light bar: Outer Flashing Lights: Switched PRIMARY / SECONDARY Center Flashing Light: Switched PRIMARY ONLY Non KKK Flashing Lights: Switched PRIM/SECONDARY	1 1 1 1

07/10/2019 13:			Page 15
PART NO 05-SY-LTMF 05-FS-0712 05-FS-0702	S	M SERIES LED WARNING LT SYSTEM - F-series LTD Flasher: None, All M Series LED's Are Programmable Flasher: Whelen ULF-44, IATS NOTE: ULF-44 Flasher for White (Alley) Lights on 5V3R Lights Only and Reds Will Flash Off of Flasher in Light Head	QTY 1 1 1
		Auto Dim Feature: Do not activate Auto-Dim Feature. NOTE: The Split M Series light will not work with the Auto Dim feature	
05-FS-10P4 11-ZZ-0909		Flash Pattern: None, Internal Light Head Flasher CAAS GVS Exception to C.8.2.12 Warning lights non-compliant to Table 1 dual mode The Warning lights pattern chosen by the Customer does not meet either the required photometric output flash or the pattern required of Table 1, Dual mode, Primary Secondary of CAAS C.8.2.12. An exception has been taken for only this element to the standard.	1 1
06-SW-PS01 11-ZZ-0909		Warning Light SWITCH: center console, Primary / Secondary CAAS GVS Exception to C.8.2.12 Warning lights non-compliant to Table 1 dual mode The Warning lights pattern chosen by the Customer does not meet either the required photometric output flash or the pattern required of Table 1, Dual mode, Primary Secondary of CAAS C.8.2.12. An exception has been taken for only this element to the standard.	1
ZE-00-0010	S	/ Flash Pattern: Action Flash 150 Specify pattern # and description for each type of light from the Whelen Web site	1
ZZ-ZZ-ZZZZ 05-LB-8005	X	/ / Opticom: Whelen 792H Opticom Unit installed Strobe Infra-red Opticom, Stand Alone unit. EMIT-792H NON waterproof system	1
05-LB-8121		Opticom Location: Built into 600 Series Center Front Warning	1
ZZ-ZZ-ZZZ		Light	1
05-PH-LSXA 05-PH-LT0E		Mirror Lights: Ford Type 1, In Black ABS Bezel Additional Lights: (2) Whelen TIRLIN3, RED LED, Chrome Hsg	1
05-PH-LA11 05-PL-LX00		Red LED Lens color for Above LED Light to be: Clear	1
ZZ-ZZ-ZZZZ		<i>I</i> -	1
05-PH-LSY5	S	/ (2) Grille Lights: Whelen 5V3R Surface Mount NOTE: Program Grille Lights #21 Action Flash	1
		NOTE: Red/White Lights to Flash in Primary Mode and Off in Secondary Mode and Puddle Lights To Activate With Parking Brake	
05-PL-LXDR		Lights: (2 Whelen 500 V Series 5V3R, LED, RED ILOS	2 2
ZZ-ZZ-ZZ01 05-PM-LE09	S	/ Location:Grille Lights / (2) Intersection Lights: Whelen 5V3R, LED, 15° Cast Housings NOTE: Red/White Lights to Flash in Primary Mode and Off in Secondary Mode and Puddle Lights To Activate With Parking Brake	1

07/10/2019 13:	:48:2	28	Page 16
PART NO	S	DESCRIPTION	QTY
05 BU I B00		Latercation Haminus (0) Oast FTIDC Linkto 45% F Carias 2047.	
05-PH-LR00		Intersection Housings: (2) Cast, 5TIR6 Lights, 15°, F-Series 2017+ HOUS-52717	1
		HOUS-52718	
		for LED TIR6	
		TOT LED TITLE	
05-PL-LXDR	S	Lighst: (2) Whelen 500 V Series 5V3R, LED, RED ILOS	2
ZZ-ZZ-ZZ 01	S	/ Location: Front Intersection Lights	2
<u> </u>		I [']	2
05-PH-LT23		DELETE (1) Center Front Warning Light	1
05-PH-LT25		Light(s) are moved into the front light bar.	1
<u> </u>		 	1
05-PH-LT24		DELETE (2) Front Warning Lights	1
05-PH-LT25		Light(s) are moved into the front light bar.	1
ZZ-ZZ-ZZZ		 	1
05-PM-LH70		(4) Upper Side Body Lights: Whelen M9 Series, LED	1
05-PH-LF14		M9 Series Flange, Chrome Plated, Std: Two Pair	1
05-PM-LRM1		Light: Whelen M9, LED, Split RED/BLUE LED/CLEAR Lens,	4
0E EC 4040		Programmable ILOS	
05-FS-1212		Opt A Programmable Lights Split Color - Color 1 Flashes Opposite Color 2	4
05-PL-LY64		LED Color for Above M9 LED : RED/BLUE	4
05-PL-L104 05-PM-LV10		LED Color for Above M9 LED . RED/BLUE Lens Color: Clear	1
05-PM-LS44		LED M9 Red - Two Pair	1
05-PM-LS44	X		1
05-PM-LI50	-	(2) Rear Intersection Lights: Whelen M7 Series, LED	1
05-PH-LF22		M7 Series Flange, Chrome Plated, Std: Each	2
05-PM-LR34		LED M7 Red/Blue Split - Pair	1
05-PM-LRC3		/ Light, Whelen LED M7 Red/Blue LED/Clear Lens - Programmable	2
05-PL-LY16		LED Color for Above M7 LED : RED/BLUE	2
05-PM-LV10		Lens Color: Clear	2 2 2 2
	X		2
05-PH-LT59		DELETE (2) Rear Warning Light	1
05-PH-LT26		Light(s) are moved into the rear light bar.	1
77-77-777			1
05-PM-LK70 05-PH-LT6B		(2) Additional Rear Body Lights: Whelen M9 Series, LED Location: REAR, (1) aligned w/ EACH upper window in RA doors.	1
05-PH-LT0B		Location: REAR, (1) aligned w/ EACH upper willdow in RA doors. Light: Whelen M9, LED, RED LED/CLEAR Lens, Non-Programmable	1 2
03-1 W-L104		IATS	
05-PH-LF10		/ M9 Series Flange, Chrome Plated, Std. Each	2
05-PL-LY54		LED Color for Above M9 LED : RED	2 2 2 1
05-PM-LV10		/ Lens Color: Clear	2
06-EL-2500		Brake Light Override: Wire into rear warning lights above	1
		LIGHTS SHALL STEADY BURN FOR BRAKES	
		No flasher included in this option code.	
			
06-EL-2502		/ Override Function: During ALL modes	1
ZZ-ZZ-ZZZZ		 DELETE (4) De en Conton Mounines Links	1
05-PH-LT69		DELETE (1) Rear Center Warning Light	1
<i>05-PH-LT26 05-PN-KPM9</i>	c	Light(s) are moved into the rear light bar. Additional Upper Side Warning Lights: (4) Whelen M9 Series, LED	1
05-PN-KPM9	3	Additional Opper Side Warning Lights: (4) Whelen M9 Series, LED M9 Series Flange Chrome Plated, Std: Pair	2 2
OG-1 11-LI 1Z		i mo deries i lange om ome i lateu, stu. Fall	

07/10/2019 13			Page 17
PART NO	S	DESCRIPTION	QTY
05-PM-LT26		Light: Whelen M9, LED, Split RED/BLUE LED/CLEAR Lens, Programmable IATS	4
05-PH-LF10		M9 Series Flange, Chrome Plated, Std. Each	4
05-PL-LY64		LED Color for Above M9 LED : RED/BLUE	4
05-PM-LV10		Lens Color: Clear	4
25-PH-LT84		LED Lights: Programmable, Single Light	4
ZZ-ZZ-ZZZ		/	2
ZZ-ZZ-ZZZ		 	1
05-ZZ-ZZZZ	C	FLECTRICAL CELL 2	1
	S	ELECTRICAL - CELL 3 NOTE: Inverter, Suction Pump, Breaker Box, Etc. Will be Mounted in M8 Compartment	
05-IL-12MW	S	/ Cab Ovhd Lights: (2) Whelen 3SRCCDCR, Red/White LED Installed NOTE: One Above Driver and One Above Passenger Seat in Cab With Switches in Cab Switch Panel	2
06-00-0005	S	Customer Supplied Part(s) : Yes See account	1
00-00-0003	3	1) AEV to Install (3) Antennas Shark Fin, Short and Longer Antenna	′
		2) AEV to Install Two Radio Heads	
		3) AEV to Install (1) Radio Trunk	
		4) AEV to Install Two Radio Cables	
07-00-0182	S	Custom design composite console included in Order !!!!!!!!!!	1
06-00-0001 <i>05-IL-2007</i>		Rechargeable Flash Lt, Fire Vulcan #44401 Orange, Installed	1
05-IL-10L0	S	Rechargeable Flash Et, Fire Vulcair #44401 Orange, installed Location: M-5 Compartment, wall No 1, NOT ON/NEAR FLOOR!	2
05-IL-10L0	X		2
05-IL-10P0		Powered: 12v Constant Hot	2 2 2 2 1
06-AL-47SU		Hand Held Spot Light: None	
06-AL-49S4		/ CAB Spot Light: GoLight #20204: LED, White Housing with Hard	1
		Wired Remote includes (1) hardwired joystick remote	
		LITE-GO20204 LED spotlight WHITE	
		No Rock guard from GO-LIGHT available for this model.	
06-AL-5000		CAB Spot Light Mount: Cab Roof	1
06-BA-3F41		Batteries: 2 - Battery System - Type I Only	1
06-BA-3F53 06-BA-3FJ1		Batteries: (2) Underhood (Type I) Battery Make: (2) OEM	1
00-BA-3F31		Battery Brand must match.	1
		Datiery Drana materialism	
06-EC-0890	S	/ Battery Switch & Auto-Battery Shut-down (Kill Circuit): 10 Minute MDUL-ss42 Non adjustable timer set for 10 minutes. This will Include Cole Hersee 2484-16 Paddle on Driver Side of Cab Console Timer Activated by ignition.	1
ZZ-ZZ-ZZZZ		I	1
06-EC-056C		POWER Door Locks: Module Doors, X-Series Units type 1 Ford	1
		2 11 2 2 2 2 2	

07/10/2019 13		Page 18
PART NO	S DESCRIPTION	QTY
06-EC-050B	Door Locks, Tied into OEM System X series Type 1 Ford	1
	Module power locks are tied to the Chassis power lock system	
	MDUL-BDLM506A Module door lock control F3/4/550	
06-EC-0512	OEM Koy Esh: Included w/Chassis Unlock Med Doors	1
00-20-0312	OEM Key Fob: Included w/Chassis Unlock Mod Doors Must also operate the Module Power Door Locks Optioned.	'
	Must also operate the Module Fower Boor Locks Optioned.	
06-EC-05B9	Door Lock Switches : Integrated in Interior Entry door Handles -x	1
06-EC-05M0	ONLY the following doors shall have power door locks:	i
06-EC-05X1	Electric Door Lock: (M-1) Left Front Compartment	1
06-EC-05X2	Electric Door Lock: (M-2) Left Middle Compartment	1
06-EC-05X4	Electric Door Lock: (M-3) Left Rear Compartment	1
06-EC-05X6	Electric Door Lock: Rear Access Doors	1
06-EC-05X7	Electric Door Lock: (M-5) Right Rear Compartment	1
06-EC-05XA	Electric Door Lock: Curbside Access Door	1
06-EC-05XB	Electric Door Lock: (M-7) Right Front Compt	1
06-EC-05XC 06-EC-081X	S Electric Door Lock: (M-8) Right Front Compt Remote Keypad #1: Installed, X-Series Units	1
00-EC-001X	MTB to cut mounting holes per provided drawing.	1
	with to dut mounting holes per provided drawing.	
06-EC-0816	Location: By Curbside Door	1
06-EC-2800	Door Unlock Switch, Momentary, Exterior, hidden	1
	NOTE: Installation of Remote Door Lock Switch feature may increase likelihood of	
	unauthorized entry into vehicle. By checking this option, purchaser further agrees	
	to hold REV or chassis manufacturer harmless for any loss of vehicle or contents	
	caused by unlawful access.	
06-EC-0802	Location: Front Grille/Bumper Area	1
ZZ-ZZ-ZZZZ		1
06-EC-3046	Cam Sys - Four White Backup, Interior bullet, Right, Left turn	1
	(1) AOM7694HD4 Kit	
	Rear camera is a fixed position White surface mount camera.	
	(1) CAMR-CAMCCD	
	(1) ADAP-31100041	
	To Consist of:	
	7" Monitor Mount Exterior White Backup Camera, fixed position.	
	Patient Area Bullet Camera	
	Right/Left Side Turn Cameras (reverse image) chrome housings	
	(2) Camera Cable (50Ft)	
	(2) Camera Cable (15Ft)	
00 50 00-0	S Monitor: Mount on Pedestal At Headliner	1
06-EC-3070	Recording System - ASA Voyager DVR DVR4100KIT Recorder	1
	(1) RCDR-DVR4100KIT	
	>> AEV/ REV installer to check SD card install after installation for recording capability of camera input and playback.	
	(1) Voyager DVR4100 KIT recorder (allows up to 4 camera inputs) also has video	
	out for displaying rear back-up events	
	(1) 32gb SD Card (Included)	
	(1) GPS antenna (Included)	

07/10/2019 13:	:48:	28	Page 19
07/10/2019 13: PART NO 06-EC-4300 06-EC-43B0 06-EC-CCX2	S	Back-up Alarm: Standard 102DB Cut Off Switch: Auto reset ,momentary style Circuit Board, 12V: Multiplex LX1 F3/4/550 X-Series STD Includes the Following Items, please add to order: All the following need to be added as usual: Rear Digital Thermostat/Control for A/C, includes front also Door Open 15, minute check out timer for 15 Oxygen Contents Gauge front and back. LOCATION NOTE: Switch Panel to be Located on Cab Console SWITCHING: Master Switch Load Manager/Sequencer Primary/Secondary Siren/Horn Left Scene Rear Scene Right Scene Rear Dome Lights Back up Alarm Disable Silent Intercom Lights Door Ajar Light Compartment Open Light Voltmeter, Digital Hands Free siren Off/On	Page 19
06-EC-1304 06-EC-CB07 06-EC-SPA7 06-EC-SPA8 06-EC-SPZ0		Howler activation switch Silent Intercom Indicators Voltage Monitor, LX-1, w/ Warning Light Circuit Protection, 12V: FET - Auto-reset LX1 Main Switch Panel Location : Patient Area LX1 Panel: SWITCHING: Dome High/Low Red Dome Fluorescent Lights Silent Intercom Switches Driver Buzzer Exhaust Fan Thermostat Control for Heat/A.C. Heat/A.C. Fan Control, 3 speed. Cabinet lighting	1 1 1 1
06-EC-SPZ1 06-EC-SPA1 06-EC-SP70 06-FS-1106 06-RR-240Z	S	Pt Area LX1 Switch Panel: One Standard In Patient Compartment Switch Panel Located: Streetside A/A Patient area switch panel mounted :flat to action wall std Sequencer, LX1, with load manager,	1 1 1 1

07/10/2019 13	:48:	28	Page 20
PART NO	S	DESCRIPTION	QTY
07-RR-020F		Add-on Console: Type 1	1
07-RR-02A0		Drink Holder: (2) In Cab Console, Big Gulp Size	2
07-RR-02A5		Drink Holder Location: Front	2
07-RR-02E0	c	Note Book Slot: Single - Full width by 3" wide	1
	S		1
07-RR-02E1	S	SPECIAL Open storage space rear of note book slot 5x5" x full	7
	_	width	
07-RR-2301	S	Customer Radio: Space for Customer installed Radio Heads	2
		radio head below fire come control	
		also Customer radio power trunk, Custoemr supplied install behind?	
		????????	
07-RR-2352		Glove Storage : Stainless Steel triple Glove box dispenser:Top	1
07-KK-2332		loader	′
		(1) triple glove box holder top load installed	
		(1) GLOV-BUTLER3	
		Brushed stainless steel	
77 77 777	0	I may not place by the car be also well of each below talls (them windows	
ZZ-ZZ-ZZZZ	S	mount glove butler on back wall of cab below talk tthru window	1
		above console.	
07-RR-24A0		Console Finish: Black, Textured "Easy Grip"	1
07- ZZ-ZZZ		/	1
06-RR-S502		Custom Wood/ Ez-Grip Console Drawing request	1
		AEV Engineering to develop drawing for console and send to Account manager for	
		approval prior to release.	
		Console to include standard items.	
		Removable top panel, inset	
		Front location for Switch panel.	
		all other items at special cost.	
07 DD 0440		L. Canada Finiah, Black Tautuwad IIFaan Cuinil	4
07-RR-24A0		Console Finish: Black, Textured "Easy Grip"	1
07- <u>ZZ-ZZZ</u>		 	1
06-EC-GR01		Ground Straps, Module to Frame: (Qty 4) Braided	1
06-IA-0001		Converter , 125V to 12V: NONE	1
06-IG-0310		Inverter : Vanner 20-1050 CUL-DC - Full Mod	1
06-EC-03AS		Portable Equip Charging Circuits: Included in Inverter	1
06-EC-03AA		PREWIRE LOCATION: (1)Cab Console, (1) Behind A/A	1
06-EC-03C9		Portable Equip Pwr Source: Ignition and/or Shoreline	1
00-20-0309			,
		Reference 3.7.7.2 and Figure 7 in KKK-A-1822E	
06-IG-03A0		Pattory Charger/Conditioner, EEA Puilt into Investor	4
		Battery Charger/Conditioner: 55A - Built into Inverter	1
06-IG-0003	_	Built-in Battery Charger: Enable - Wire to Batteries	1
06-IG-04L0	S	,	1
06-IG-04M0		On Floor near wall #1	1
06-IG-04C0		Lexan Cover: Over Inverter	1
ZZ-ZZ-ZZZZ		· 	1
06-MC-0901		Low Voltage Buzzer: Installed, Multiplex units	1
		COMMUNICATION RADIO(S) RELATED	
06-00-0001			'
00-00-0001			
00 50 0040		RADIO POWER]]
06-EC-03A0		Radio Power No 1: 30A, Pos and Neg, 10 awg Wires	1
06-EC-03B0		Radio Power Source: Battery Switch Hot	1
06-EC-03J0		LOCATION: Inside Cab Center Console	1
06-EC-03A2		Radio Power No 2: 30A, Pos and Neg, 10 awg Wires	1

07/10/2019 13	:48:29	Page 21
PART NO	S DESCRIPTION	QTY
06-EC-03B0	Radio Power Source: Battery Switch Hot	1
06-EC-03E0	LOCATION: Behind Action Area Board	1
06-EC-03A3	Radio Power No 3: 30A, Pos and Neg, 10 awg Wires	1
06-EC-03B0	/ Radio Power Source: Battery Switch Hot	1
06-EC-03E0	/ LOCATION: Behind Action Area Board	1
06-EC-03A4	Radio Power No 4: 30A, Pos and Neg, 10 awg Wires	1
06-EC-03B0	/ Radio Power Source: Battery Switch Hot	1
06-EC-03J0	/ LOCATION: Inside Cab Center Console	1
06-EC-03A5	S / Knox Box Power Prewire: 30A, Pos and Neg, 10 awg Wires	1
06-EC-03B0	/ Radio Power Source: Battery Switch Hot	1
06-EC-03F0	S / LOCATION: M-8 Compartment	1
ZZ-ZZ-ZZZZ	, 	1
	ANTENNA LEADS	1
06-RR-0100	S GPS Transmitter Antenna Antenna, No 1: Customer Supplied	1
	Shark Fin Shaped Antenna	
	·	
06-RR-01O1	ORIGINATION POINT: Roof Port No 1	1
	Roof Port No 1 is lined up with the edge of the curbside entry door, and centered	
	side to side as possible	
	·	
06-RR-01P0	Port Plate: None - Access thru center upholstery panel	1
06-RR-01T7	TERMINATION POINT: Center Console w/ 36" Tail	1
06-RR-1712	/ Antenna Base w/ Coaxial Cable: KE794 #2	1
	MOUN-KE794	
06-RR-0102	ORIGINATION POINT: Roof Port No 2	1
	Roof Port No 2 is located approximately 24 inches back rearward of roof port #1.	
06-RR-01P0	Port Plate: None - Access thru center upholstery panel	1
06-RR-1772	Termination Point: Inside Center Console	1
	MOÚN-K794	
06-RR-1712	X Short Antenna, No 2: Customer Supplied	1
	NOTE: Install Shorter Antenna on KE794 Base to be Long Enough to Reach	
	to Cab Console. Customer Supplied Coax Cable Was not Long Enough.	
06-RR-1714	Antenna Base w/ Coaxial Cable: KE794 #3	1
	MOUN-KE794	
06-RR-01O3	ORIGINATION POINT: Roof Port No 3	1
	Roof Port No 3 is located approximately 24 inches back rearward of roof port #2.	
06-RR-01P0	Port Plate: None - Access thru center upholstery panel	1
06-RR-1772	Termination Point: Inside Center Console	1
	MOUN-K794	

07/10/2019 13	48:29		Page 22
PART NO	S	DESCRIPTION	QTY
06-RR-1714	NOTE	Long Antenna, No 3: Customer Supplied It Install Longer Antenna on KE794 Base to be Long Enough to Reach To Console. Customer Supplied Coax Cable Was not Long Enough.	1
06-RR-1100	NOTÉ Behin	traction - (2) Radio Cables, Customer Supplied prior to production, Installed : Install One Longer Radio Cable From Radio Head in Action Area to d Passenger Seat to Radio Trunk and Install One Shorter Cable From Head in Console to Radio Trunk Behind Passenger Seat.	1
06-RR-1809 <i>06-RR-1910</i>	Install Cab C	Patient Area Radio: None Intercom, Kwik Kit, Havis, Push To Talk system Constant Hot Box in Action Area and Install Push to talk box in Center of seiling NTE-PT	1 1
ZZ-ZZ-ZZZZ 06-SO-0000	<i>J</i> -	125V SHORE LINE AND OUTLETS	1
06-00-0001 06-SO-0500 06-SO-00L1	S /-	Shore Line Inlet: 20A Super Auto Eject, ILOS Inlet location: Left Rear of Modular Box - In Kick Plate Drawing"	1 1
06-SO-0805		tor light located above shore line inlet.	1
06-SO-1005 ZZ-ZZ-ZZZZ 06-SO-10TT	/- NOTE	Cover, Yellow, Shore Line Inlet: 20A Super Auto Eject, STD Inpower Timer: VCM-05-01SF, Installed This allows the auto eject to be wired to the ignition switch In lieu of g into the OEM starter circuit and is required any this vehicle	1 1 1
06-SO-0901	Requi	Shore Line Inlet: 30A, Auto Eject, IATS res additional space behind interior cabinets for installation. Install Shore peration label.	1
06-SO-00L1 06-SO-0805	-	Inlet location: Left Rear of Modular Box - In Kick Plate Indicator Light, Shore line: Weather proof, power by SL, Red tor light located above shore line inlet.	1
06-SO-1080 ZZ-ZZ-ZZZZ ZZ-ZZ-ZZZZ 06-SO-1400 06-SO-1401	 	**125 Volt OUTLETS** - 125 VAC Outlet, No 1: 15A, Hospital Grade, IVORY 5 VAC outlets shall be back lighted when power is applied to the outlet.	1 1 1 1
06-SO-14L1		LOCATION: Action Area, standard location	1

07/10/2019 13			Page 23
PART NO	S	DESCRIPTION	QTY
06-SO-14O2		Outlet mounting ORIENTATION: Vertical	1
ZZ-ZZ-ZZZZ			1
06-SO-1402		125 VAC Outlet, No 2: 15A, Hospital Grade, IVORY	1
06-SO-14L1		LOCATION: Action Area, standard location	1
06-SO-14D2		Outlet mounting ORIENTATION: Vertical	1 1
ZZ-ZZ-ZZZZ		Outlet mounting ONIENTATION. Vertical	
06-SO-1403		125 VAC Outlet, No 3: 15A, Hospital Grade, IVORY	1
06-SO-14L3		LOCATION: RF ALS, (See Drawing)	1
06-SO-14O2		Outlet mounting ORIENTATION: Vertical	1
<i>ZZ-ZZ-ZZZZ</i>		/	1
06-SO-1404		125 VAC Outlet, No 4: 15A, Hospital Grade, IVORY	1
06-SO-14L4	S	LOCATION: Special Location - (See Drawing)	1
		LOCATION NOTE: Foot of squad bench, forward facing, near portable	
		suction.	
06-SO-14O2		Outlet mounting ORIENTATION: Vertical	1
ZZ-ZZ-ZZZZ		•	1
06-SO-1100		/ **INTERIOR 12 Volt OUTLETS**	'
06-SO-10XZ		12 Volt outlets two separate single outlets	1
06-SO-10Y0		12V Outlet, No 1: Power Point Single - Wire thru Med Isolator	1
06-SO-11L1		LOCATION: Action Area, standard location	1
<i>06-SO-1402</i>		Outlet mounting ORIENTATION: Vertical	1
06-SO-1910		Power Source: Medical Isolator , Batt Sw Hot	1
<i>ZZ-ZZ-ZZZZ</i>		/	1
06-SO-1102		12V Outlet, No 2: Power Point - Wire thru Med Isolator	1
06-SO-11L3		LOCATION: RF ALS, (See Drawing)	1
06-SO-14O2		Outlet mounting ORIENTATION: Vertical	1
06-SO-1913		/ Power Source: The SAME as outlet No 1	1
ZZ-ZZ-ZZZ		/ 	1
06-SO-1103		12V Outlet, No 3: Power Point - Wire thru Med Isolator	1
06-SO-11L8		LOCATION: In Cab Console	1
06-SO-14O2		Outlet mounting ORIENTATION: Vertical	1
06-SO-1913		Outlet mounting ONENTATION: Vertical Power Source: The SAME as outlet No 1	1
ZZ-ZZ-ZZZZ		Fower Source. The SAME as outlet NO T	1
	0		1
06-SO-116K	3	12V DUAL USB 4.8A,. (3) Outlets 2.4 A, switch panel large hole mount	3
		Mounts in Oversize switch panel hole,	
		Includes pilot light to indicate power to connection.	
		Open slots, no dust or water closures.	
		4.8 Amp max output total.	
		DUAL USB outlet	
26-SO-11L5	S	LOCATION: (1) In Cab Console, (1) In A/A & (1) On C/S Wall	3
ZZ-ZZ-ZZZ		 	3 3 1
06-SS-070A		Siren: Whelen, WS295HFSC9, F-Series	
06-SS-SA27		Siren Speakers:CPI: SA4327 Through the bumper 2017 F series STD	1
06-SS-SW01		/ Siren / Horn Switch: In Cab Console	1
		This standard rocker switch in the front switch panel shall select horn ring output:	
		Siren Functions OR OEM Horn	
06-SS-SW04		microphone clip not installed Tie to Siren Microphone cord.	1
06-SS-1313	2.	Air Horn System: BUELL - F1	1
20 30 .010	9	NOTE: The foot air horn switches shall be interfaced with the door light	'
		switch to prevent air horn from being activated when occupants are	
L		onition to provent an norm from some activated when occupants are	

	29	Page 24
S	DESCRIPTION entering/dismounting vehicle.	QTY
S		1
	NOTE: The foot air horn switches shall be interfaced with the door light switch to prevent air horn from being activated when occupants are entering/dismounting vehicle.	
	Air Tank Mounted to Frame Rail Additional Horn Discharge ENABLE: Both Cab Doors Closed Horns will NOT discharge with either cab door open.	1
		1 1 1 1 1
	Siren, Whelen Howler, IATS, w/Bracket - F4/550 2017+ > Not available for F-350 4x2 due to ground clearance issues Must mount off center, not in the center air flow area.	1
X	Howler, Mounting Bracket: F3/4/550 - HWLRB22 - 2017 Siren Pre Wire for Q2B Mechanical Siren NOTE: Install Pre Wire in Front Bumper Area for Dealer to Install Q2B Siren with Solenoid and Rocker Switch in Console. Make sure that prewire has equal size wires	1 1
•	CABINET SHOP - Modular	1
5	j	1
	ALL DIMENSIONS ON ORDER WILL HAVE A VARIABLE OF + OR - 1/4" Unless specified as a minimum dimension.	1
		1
	Rubber Cactus Mat : On Cabinet Floors, Drawers & Shelves, X-Series	1 1
	units Aluminum Cabinet Backer-Sound Barrier # INSU-KSU on back side of aluminum cabinets	1
S	Interior Color Packages Mica Colors: Gloss Gray over Stainless Steel Parting Line: C/S: Top of squad bench including rear filler. S/S: Bottom of Cabinet F Stainless Steel on above areas only. No other cabinet gets stainless, unless specified otherwise below.	1 1
	s s s	entering/dismounting vehicle. S

07/10/2019 13:	:48:	29	Page 25
PART NO	S	DESCRIPTION	QTY
07-AS-1200		Stainless Steel - FACE of Squad Bench including rear filler	1
		End of squad bench shall be MICA.	
07.40.4400			
07-AS-1400		Stainless Steel - Lower Face of Wall Cabinet - Included	1
ZZ-ZZ-ZZZ		 Debugget and to Toma (Onlaw Louis Col FAR Colours Lottels	1
07-00-PC0N		Polycarbonate Type/Color: Lexan - CLEAR Secure Latch	1
07-00-PH01	C	Window Handles: Full Length Extruded	1
07-00-PS02 07-00-PS02	S X	Security Tag Holes: Custom With Secure Latch Windows	1
11-SX-0400	^	 Attendant Seat EVS Child 6pt 6degree with Black Belt X series ILOS	1
11-3X-0400		Attendant Seat Evs Crind opt odegree with Black Belt X Series iLOS	'
11-SX-0401		Attendant Seat:EVS, 6 deg Child Safety Seat, 6-Pt, Gun Metal, x	1
11-02-0401		series ILOS	′
		ARMRESTS are not permitted on this seat as of 02-15-16	
		SEAT-14VEV6G6T	
		SEAT-LOGO AEV REV	
07-SE-0901		BASE: Swivel base for EVS Child Seat Option	1
		(2) Position Swivel, Forward or Rearward, Only	
		BASE-1882S	
		T-Handle Standard at AEV 10/01/18	
07-AC-4500	S	AC CABINET: Evaporator,Location : Above RF ALS X Series STD	1
		composite	
		Composite Construction	
07 DV 0000		D	
07-BK-0002	S	Bulkhead Cabinet Above Talk Thru: Installed IATS, See Drawing	7
		NOTE: Cabinet Above Talk Thru With Medix Safe Installed on the Pass Side of Cabinet "J"	
		SEE DRAWING	
		SEE DRAWING	
07-DR-LXX8		Doors, Dual, Overlay Hinged, 3/8" Lexan, Small Cabinet	1
07-HW-SX11		Round Pull Latch: Non-locking - Chrome Finish	1
11-X0-0010		/ Stowage rating label - Southco round latch 10 pounds applied	1
		each	
07-ZZ-ZZZZ			1
07-CA-03XA		CURBSIDE UPPER: Over Squad Bench	1
07-CU-X107		Cabinet K1 & K2: (2) Openings with flush center divider	1
07-DR-LXX2		(K) Door: Single Flip Up 3/8" Lexan	2
07-DR-LZ20		Hinges: Reel Torque	2
07-HW-SX11		Round Pull Latch: Non-locking - Chrome Finish	2
11-X0-0010		/ Stowage rating label - Southco round latch 10 pounds applied	2
07.74.0000		each	
07-XA-3800		Divider System For Adjustable dividers w/ ABS slotted track -	2
07 VA 2040	C	Each Cabinet	0
07-XA-38A0	S	/ Lexan Dividers, (4) In Each Cabinet NOTE: Four Lexan Dividers in Cabinet K1 and Four Lexan Dividers in	8
		Cabinet K2.	
		Cabinet NZ.	
07-ZZ-ZZZZ		I	1
07-CA-04A8		/ / Cabinet, (3) Glove Storage, over C/S Entry Door	1
07-DR-LX04		Cabinet, (3) Glove Storage, over 6/3 Entry bool Door: Overlay He, 3/8" Lexan - 3-glove dispense thru	1
OF DICEROT		, 2001. Overlay 110, 6/6 Lexail o glove disperse dila	

07/10/2019 13:	:48:	29	Page 26
PART NO	S	DESCRIPTION	QTY
07-HW-SO11		Round Pull Latch: Non-locking - Chrome Finish	1
11-X0-0010		Stowage rating label - Southco round latch 10 pounds applied	4
11-20-0010			,
		each	
07-ZZ-ZZZZ		<i> </i>	1
07-CA-0400	S	CURBSIDE REAR: Rear of S/B, Mods	1
01-0A-0400	0		•
		NOTE: Cabinet to be Located Above M5 Compartment and Above Drawers	
		and Will Have Dual Lexan Doors to Open Towards Aisle.	
07-CU-K103		Cabinet L : (1) Opening with (1) Adj shelf, shelf std.	1
	_		
07-DR-LX24	S	Doors: Dual Overaly 3/8" Lexan (1) Hinged Right and (1) Hinged	7
		Left	
07-HW-SO11		Round Pull Latch: Non-locking - Chrome Finish	2
11-X0-0010		•	
11-20-0010		Stowage rating label - Southco round latch 10 pounds applied	
		each	
07- ZZ-ZZZ		 	1
07-ZZ-ZZZZ		, 	1
-		L DE ALC Cohinete Aluminum Comments Comments	
07-RF-4500		RF ALS Cabinet: Aluminum Composite Construction	1
		X-Series Design with 45°	
07-CA-04AY		X series Angled Glove at RF ALS not required	1
			<i>!</i>
07-RF-0003		Interior Color: White Anti-Microbial Powder Coated	1
07-RF-I1X0	S	Cabinet I-1: Deleted for Installation of HVAC system	1
07-ZZ-ZZZZ		Í	1
07-RF-I2X0		Cobinet L2: Standard Aluminum Composite Construction	ا ا
		Cabinet I-2: Standard, Aluminum Composite Construction	!
07-DR-LXXC		(I) Doors, Dual, Solid - RF ALS Cabinet	1
07-HW-SO11		Round Pull Latch: Non-locking - Chrome Finish	2
11-X0-0010		/ Stowage rating label - Southco round latch 10 pounds applied	2
11-20-0010			_
		each	
07-RA-IOA1		Outside Access: Thru M-7 (RF) Compartment door.	1
07-RF-0001		Interior MICA Color: White cabinet liner	1
07-TS-12X2		Shelf Track: Small Alum. Unistrut, X-Series, RF ALS cabinet	1
07-CA-06VD			1
		(2) Shelves: Adjustable, Aluminum, White Anti-Microbial coated, STD	
ZZ-ZZ-ZZZZ			1
07-RF-I30D		Cabinet I-3: DELETE, Combine with I2 Cabinet	1
07-ZZ-ZZZZ		, 	1
07-ZZ-ZZZZ		1	4
		T	
07-RR-0064		Right Rear Cabinet: Cover over M-5 Compartment	1
07-ZZ-ZZZZ			1
07-SA-0000			1
		L COLLAD DENICLE Alternia um Comercacita Constitucation V Conic	
07-SB-4500		SQUAD BENCH: Aluminum Composite Construction X Series	1
07-CA-42A0	S	Tilt out - See Drawing	1
		tilt out door at squad bench	
		Trim opening	
		non locking round latch to secure	
07-HW-SO11		Round Pull Latch: Non-locking - Chrome Finish	1
		•	
11-X0-0010		/ Stowage rating label - Southco round latch 10 pounds applied	1
		each	
07-SB-1003	S	LOCATION: At head of Squad Bench	1
	-	NOTE: Tilt-Out Towards Aisle.	
		ear lettalae/herei	
07 CD 4004	_	L CHARRE CONTAINED, NOVE	
07-SB-1004	S	SHARPS CONTAINER: NONE	1
0. 02 .00.			

07/10/2019 13:	<u>48:</u> 2	29	Page 27
PART NO	S	DESCRIPTION	QTY
07-DW-0010	S	Drawers: (2) Aluminum Facing Front of Truck Above M5	1
		Compartment	-
07 DD WD45		•	2
07-DR-WD15		Drawer Front: Flush Fitted	2
07-CA-2100		TRIM: U-shaped Door, J-trim opening	2
07-HW-SO11		Round Pull Latch: Non-locking - Chrome Finish	2
11-X0-0010		/ Stowage rating label - Southco round latch 10 pounds applied	2 2
		each	
07-HW-SL01		Drawer Slides: Full Extension, Detachable 130 LB capacity STD	2
07-SB-100N		/ SHARPS CONTAINER: BD 5.4 Quart Sharps Wall mount with	1
OF OB TOOM		locking cradle IATS	,
		CONT-305517	
		BRKT-305447	
07-SB-0996	S	LOCATION: On C/S Wall - See Photo on file	1
07-SB-0001		See Upholstery Section for Post and Wheel Cup info	1
07-SB-1001	S	Storage Under Lid - Configure to M-5 Compartment Size	1
07-SB-1001	9	Storage Order Eta - Comigare to in-3 Compartment Size	1 1
07-SB-LH00		Hinge, Squad Bench Lid(s): Butt Style Hinges	1
07-SB-LH03		Lid Checks: Gas shock, Dual Action	2
07-SB-LH07		Latch, Squad Bench Lid: Slam Action Paddle, W keeper Compliant	2
		J3058	
11-X0-0030		Stowage rating label - Squad bench interior, entire area 80 lbs	2
07-TR-SB01		Edge Trim, Lids: Band w/ Laminate and J-Trim Protection	1
ZZ-ZZ-ZZZZ		Lago Timi, Elao. Bana w/ Laminato ana o Timi i Totootion	1
		Postroint Not Detectable at boad of C/D Valley Wathing II OC	
08-MH-1401		Restraint Net, Detachable, at head of S/B, Yellow Webbing ILOS	7
		Seven (7) NET Attachment Points:	
		Yellow in color.	
		NETS-MOD2Y	
08-ZZ-ZZZ			1
07-TC-4500		TOP CABINETS, - w/ 72" HR, Aluminum Composite Construction X STD	1
07-TC-A450	S	Cabinet A: Delete due to Duo Therm	1
	3	Cabinet A. Delete due to Duo Therm	1 1
07-ZZ-ZZZZ			
07-TC-B450	S	Cabinet B: Non angled Extended Over A/A (72"HR) X Std	1
07-CB-DR75		Doors; Cabinet B Secure Latch Restocking Sliding Window Top	1
		hinged x	
11-X0-0016		Stowage rating label Secure Latch Restocking Sliding Window 30	1
		pounds applied	
07-CA-06VA	S	(2) Shelves: Fixed, Aluminum, White Anti-Microbial coated	1
07-CA-00VA	3	One on Each Side of Divider	<i>'</i>
		One on Lacii Side of Divider	
07 77 777			
07-ZZ-ZZZZ		_ 	1
07-WC-0000			1
07-WC-4500	S	WALL CABINET: No CPR Seat w/Extended A/A, Aluminum	1
		Composite Construction X	
		NOTE: Extend Action Area Rearward Approx. the Length to the Middle of	
		the Standard Placement of the CPR Seat Location.	
		and diamaka i kacement of the of it ocat Location.	
07 14/0 4 43/4	_	L. Antion Avon Entondard	
07-WC-AAX1	5	Action Area: Extended	1
07-CA-480V		1/2" Solid Surface A/A Tray No Bio Waste, X-Series ILOS	1
		No sharps or trash cutouts requested.	

07/10/2019 13:	<u>48∷</u>	30	Page 28
PART NO	S	DESCRIPTION	QTY
07-CA-483Z		Gibraltar (Solid Surface) Material COLOR: Zen Grey 9115GS	1
07-CA-4U08		Solid Surface extension in front of O1 X series std	1
07-DW-0004	S	Drawer, Aluminum Report Desk: Report supply storage under	1
07 277 0004	0	desk	′
07-DR-WD14			1
		Drawer Front: Flush Fitted	1
07-CA-2100		TRIM: U-shaped Door, J-trim opening	1
07-HW-SO11		Round Pull Latch: Non-locking - Chrome Finish	1
11-X0-0010		Stowage rating label - Southco round latch 10 pounds	1
		applied each	
07-HW-SL01		Drawer Slides: Full Extension, Detachable 130 LB capacity STD	1
07-WC-AABD		Action Area Board: Black Grip See Station No 6 for content	1
07-ZZ-ZZZ		,	1
07-ZZ-ZZZZ			1
07-WC-O10X		Cabinet O1: Standard, Aluminum Composite Construction	1 1
07-DR-LXX5		(O1)Door: Single Overlay Lexan Hinged Right	
07-DR-LXX9		Hinges: Stainless Spring Loaded	
		, ,	
07-HW-SX11		Round Pull Latch: Non-locking - Chrome Finish	
11-X0-0010		Stowage rating label - Southco round latch 10 pounds applied each	1
07-TS-12X0		Shelf Track: Small Alum. Unistrut, X-Series	1
07-CA-06VB		(1) Shelf: Adjustable, Aluminum, White Anti-Microbial coated	1
ZZ-ZZ-ZZZZ			1
07-WC-C64X	S	Cabinet C: DELETE due to increased height M-3	1
07-ZZ-ZZZZ		·	1
07-WC-D64X	S	Cabinet D1 & D2: Non Angled, Alum Composite Construct	1
07-CB-DR7D	S	Doors; Cabinet D1 Secure Latch Restocking Sliding Window	1
0. 02 22		Top hinged x	
		(1) WIND-CFTG	
		(1) WIND-OI 10	
07 TC 40V0	0	1 (4) Chalf Fixed Alimpiana White Anti Microbial Control	
07-TS-12X0	S	(1) Shelf, Fixed Aluminum, White Anti-Microbial Coated	1
11-X0-0016		Stowage rating label Secure Latch Restocking Sliding Window 30	1
	_	pounds applied	
07-CB-DR7F	S	Doors; Cabinet D2 Secure Latch Restocking Sliding Window	1
		Top hinged x	
		(1) WIND-CFTG	
11-X0-0016		/ Stowage rating label Secure Latch Restocking Sliding	1
		Window 30 pounds applied	
07-TS-12X0	S	(1) Shelf, Fixed Aluminum, White Anti-Microbial Coated	1
07-ZZ-ZZZZ	_		1 1
07-WC-E0AB	9	/ Cabinet E, (3) Glove Vertical Storage, ILOS	1
07-WO-LOAD	0	Requires Dogleg in M-3. AEV to provide MTB with engineering drawing.	'
		Requires bogieg in M-3. ALV to provide MTB with engineering drawing.	
07 DD 1 V06	_	L Doore Overley hings 2/011 Leven 2 glave dispense thru	4
07-DR-LX06	S	Door: Overlay hinge, 3/8" Lexan - 3-glove dispense thru	1
07-HW-SX11		Round Pull Latch: Non-locking - Chrome Finish	1
11-X0-0010		Stowage rating label - Southco round latch 10 pounds	1
		applied each	
07- ZZ-ZZZ		/	1
07-WC-F64X		Cabinet F: Standard, Aluminum Composite Construction	1
07-CB-DR7F		Doors; Cabinet F Secure Latch Restocking Sliding Window Bottom	1
		hinged x	
		Bottom hinged, Damper assist.	
		(1) WIND-CFBD custom sized	
		(1) WITED OF CUSTOIT SIZED	
1			1

07/10/2019 13		age 29
PART NO	S DESCRIPTION	QTY
07-TS-12X0	S (1) Shelf, Fixed Aluminum, White Anti-Microbial Coated	1
07-XA-3800	S Divider System For Adjustable dividers w/ ABS slotted track - DELETE	1
07-XA-38A0		8
U7-XA-30AU	S / Lexan Dividers, (8) Total DELETE NOTE: Four Lexan Dividers Above Fixed Shelf and Four Lexan Dividers	0
	Below Fixed Shelf	
	Bolow Fixed Gilon	
11-X0-0016	Stowage rating label Secure Latch Restocking Sliding Window 30	1
	pounds applied	•
07-ZZ-ZZZZ		1
07-WC-SS09	/ CPR Side Seat: None, See Work Order for Specifics	1
ZZ-ZZ-ZZZ		1
08-SE-1702	P4 - 4-Point PER4MAX Restraint System - with replacement indicator	1
	W/CPR	_
08-SE-1820	XS (2) on Squad Bench, (0) CPR Side Seat Per4Max Black Belts	1
	Belt retractors must mount below the top of the seat.	_
	BELT-F151390	
08-SE-SB70	S/B: (3) Sec patient restraints - 9" Sleeves Face of Bench	1
08-SE-SB76	S/B: (3) Sec patient restraints - 5" Sleeves By Hinge	1
00-20-0072	Per4Max Four-Point Seatbelts: Accept and Install Full Feature	1
ZZ-ZZ-ZZZZ		1
ZZ-ZZ-ZZZZ	[']	1
07-ZZ-ZZZZ		1
	CABINET SET - CELL 4 - MOD	1
08-00-00CX	Sub Floor, 3/4" Composite, STD	1
ZZ-ZZ-ZZZZ		1
08-00-FL32	Flooring: Lonplate II - #421 Mica (Light Gray)	1
08-00-FL01	Flooring Main Edge: 3" Recessed (1/2" deep) roll-up	1
08-AS-0005	Rear Threshold, Covered with Yellow Speedliner 6" W x Full W at rear	1
	doors	
08-AS-0006	C/S Stepwell Threshold, Polished Diamond Plate	1
ZZ-ZZ-ZZZ	/	1
08-ZZ-ZZZ	<i>J</i> '	1
08-CS-0001	COT MOUNT HARDWARE: (Full Size Mod)	1
08-CR-09PL	/ Cot Mount: Stryker, No 6390, Power Load, Installed	1
	Not Compatible with cots other than Stryker!	
	This item priced for Domestic USA usage only.	
	>>> Caution, Power load head consumes floor space toward attendant seat.	
	Spacing review between seat and cot is recommended. >>	
08-CR-09SM	No Bolt On Stryker Antler Required.	1
08-CR-1002	12v power feed location FRONT for Stryker power load 6390	1
	The 12 volt feed will be drilled and fed through the template FRONT hole for a	
	Stryker power load system. The 12volt power shall be brought through the floor,	
	with Heat shrink connectors installed.	
08-CS-00PM	•	1
	Maintain 12" between Cot and S/B.	
		_
08-CS-2000	Primary Cot position under floor reinforcement	1
	Install 8 inch wide by 5 long 1/2" aluminum reinforcement at primary cot position.	
	(Primary cot position only - REV Std)	

07/10/2019 13	:48:	30	Page 30
PART NO	S	DESCRIPTION	QTY
		**If custom location is desired, must be specified below.	
08-CS-3009		Cot mount set up for: Stryker Power-Pro 6506	1
08-CS-19A5		Cot Stop, Block: Stryker Power/Performance load. Install hook with	1
		kit	
ZZ-ZZ-ZZZZ			1
08-CS-19U1		Q Straint System: NONE	1
08-OS-04SX 08-OS-0401		OXYGEN / AIR / VACUUM System: X-Series / Oxygen Outlet No 1: Amico Console - Ohmeda/Ohio Diamond Style	
00-03-0401		Part No: OUTL-AOB	'
		Tarrivo. Cote Aob	
08-OS-040A		LOCATION: Action Area	1
08-ZZ-ZZZ			1
08-OS-0402		Oxygen Outlet No 2: Amico Console - Ohmeda/Ohio Diamond Style	1
08-OS-040A		LOCATION: Action Area	1
08-ZZ-ZZZ		j	1
08-OS-0403		/ Oxygen Outlet No 3: Amico Console - Ohmeda/Ohio Diamond Style	1
08-OS-040D		LOCATION: Curbside Wall, over the head of the S/B	1
08-ZZ-ZZZZ		/	1
08-OS-0404		Oxygen Outlet No 4: Amico Console - Ohmeda/Ohio Diamond Style	1
08-OS-040D		LOCATION: Curbside Wall, over the head of the S/B	1
08-ZZ-ZZZZ			1
08-OS-0405	_	Oxygen Outlet No 5: Amico Console - Ohmeda/Ohio Diamond Style	1
08-OS-040C	S	LOCATION: Street Side Ceiling Panel, Over Head/Chest	1
08- ZZ-ZZZ		!	1
ZZ-ZZ-ZZZZ ZZ-ZZ-ZZZZ		 	
08-OS-090X		ا Port Cyl Rack No 1: Zico QR-D-2 - D or Jumbo D - No Strap Req	
00 00 000		Through bolt the rack to metal.	'
		Through bolt the rack to metall	
08-OS-0905		Port Cyl Rack No 2: Zico QR-D-2 - D or Jumbo D - No Strap Req	1
		Through bolt the rack to metal.	
ZZ-ZZ-ZZ 01	S	Location: (2) In Curbside Step Well	1
ZZ-ZZ-ZZZZ		⁻	1
08-OS-1102		Rack No 1: Zico #QR-MV, ILOS Compliant	1
08-OS-110B	S		1
		Installed in M1Compt, Transverse	
		(1) BRKT-QRSOFD only, no oxygen restraint bracket included.	
		DELETE LATCH HOLE FOR SLIDE OUT RACK WHEN RACK IS IN THE OUT POSITION	
		POSITION	
08-OS-3605	S	O2 Regulator Mount : 5' Stainless Steel Braided Line	1
30-00-3003	3	#REGU-PTAIL	
		Add Braided line and tie it into regulator and O2 line for last section in	
		compartment.	
		•	
08-OS-1323	S	Rack Location: M-1(LF), Transverse	1
08-OS-132M		/ Set up For M cylinder ´	1
08-OS-1481		Zico oxygen tank restraint pull to tighten restraint straps (4)	1
		included	
		included with purchase of QR-MV tank holderas standard November 2017, as	
		SAE J3043 compliant.	

07/10/2019 13			Page 31
PART NO 08-OS-1712	S	I Regulator, Oxygen: Has LX1 Transducer, Inspected and Installed Installed on O2 line for pressure test, Regulator and Transducer connected and left installed on Oxygen hose. Tagged with Inspector signature/Date Protected for shipment and left on the hose, Plastic restraint tie to the oxygen bracket for shipment > Required under Change notice 11, July 1, 2018.	QTY 1
08-OS-1902 08-OS-1902 08-OS-1903 ZZ-ZZ-ZZZZ 08-OS-2300	X X S		1 1 1 1 2
08-OS-230R		<i>J Collision Restraint Strap: ZICO No 1010-120-000</i> STRP-1010	2
08-OS-23L1 08-OS-23L1 ZZ-ZZ-ZZZZ 08-OS-35ST 08-OS-3501	s x	Location: Bolt (1) to Corner of Wall #1 and Wall #2 of M3 Compartment Location: Bolt (1) to Corner of Wall #2 and Wall #3 of M5 Compartment Vacuum System: SSCOR regulator/gauge panel in A/A Collection Canister w Clip: Bemis, 1200 CC Capacity -J3043 retention compliant LOCATION: In Action Area, to the left of A/A console box. with SSCOR canister retention clip	1 1 2 1 1
08-OS-3502 08-OS-3503 08-OS-35L4 ZZ-ZZ-ZZZZ 08-OS-35VX	s s	/ [']	1 1 1 1 1
ZZ-ZZ-ZZZZ 08-ZZ-ZZZZ 09-00-001A 09-00-002S 09-00-0001 09-00-000Z 02-BT-C85V	S	TRIM / DOOR ASSEMBLY - CELL 5 Soft Touch Trim Kit: None Interior window and aluminum trim with standard satin finish- std ********Door Department X Series > New Handles, Ext: Tri-mark 1875, Free Float, Polished CNNC finish X-Series Key M-2.5 Differently	1 1 1 1 1 1 1

07/10/2019 13:		Page 32
PART NO 02-BT-C85A	S DESCRIPTION Interior Release: All Entry Doors, with bezel Emergency Access All three Entry Doors, top and bottom latches.	QTY 1
09-B1-ED0X 08-MH-0020 08-MH-0318	Entry Door Panels, Windows and Hardware, X-Series Units Interior Grab Handle Color: Yellow Antimicrobial Grab Handle, CS Entry: (1) 3-pt "L" Shaped 90, Yellow Antimicrobial New 5-12-16 Install CS door handle after cabinet installation. Close the doors and check for pinch points or cabinet interference.	1 1 1
08-MH-03E8	Grab Handles, Rear Access: (2) 3-pt "L" Shaped 90, Yellow Antimicrobial	1
09-AS-2402 09-AS-2535 09-E1-23B1	Door Panels: Mica / Upholstery / Smooth Aluminum Curbside Lower Door Panel: Smooth Aluminum Multi-Color Printed LOGO: (3) Entry Doors, Customer Logo, requires clear image Printed Multi-color logo onto panel, Back side lighted with LED lighting. ENG > to VENDOR	1 1 1
09-WI-02BX 09-WI-06BX 09-WI-09BX	"AUTO SHADE" Windows: Upper R/A Doors, Fixed Glass, X-Series STD "AUTO SHADE" Window: Upper C/S Access, Fixed Glass, X-Series STD Inverter, "Auto Shade" Windows: Part of X series windows part MDUL-WINDOW	
09-IN-3STX	Insulation: All Module Doors, 2" Black Foam, X-Series units #INSU-ARM2	1
09-IN-3STY	Insulation: All Module Doors, Black Sound Proofing, X-Series units #INSU-KSU	1
09-MH-2531	/ Grab Rail, (1), 18" Yellow Antimicrobial Rear Entry Assist std. Located: Module interior at the rear entry, mounted to the M5 compartment filler. Rail is to be through bolted into M5.	1
09-DR-WT00 09-DR-WT20 09-FM-10X3 09-FM-11X0 ZZ-ZZ-ZZZZ 07-ZZ-ZZZZ 09-ZZ-ZZZZ	Talk Through Window: Sliding Lexan Window - CLEAR Sliding Window Locking Pin: 1/4" with Lanyard Dri-Dek Compt Floors and shelves, Black, STD X-Series Dri-Dek EDGES: Black, Beveled frontal edges	1 1 1 1 1 1
	ELECTRICAL Point to Point - CELL 6	1
06-EC-1390 06-EC-1419	IV Warmer No 1: NONE IV Warmer No2 None required	1
10-00-0100	S Customer Radio: Customer Supplied Radio Head Installed in Action	1
10-MC-015F	Area S / Narcotics Box: Medixsafe Model "M2" With Wifi : Installed Dimensions: 10" H x 8"W x 12"D	1

07/10/2019 13:	48:	31	Page 33
PART NO	S	DESCRIPTION	QTY
		NOTE: Wired 12V Constant Hot and Data/Download Port to be Accessible	
		Without Having to Remove Safe. Plug Cable in Back of Safe and Run Cable	
		Out Top Left Corner of Cabinet J with Grommet.	
		Out 10p Left Corner of Cabinet 5 with Grommet.	
<i>ZZ-ZZ-ZZ</i> 01	S	Location: Cabinet Above Talk Thru Window, On Pass Side of	1
		Cabinet "J"	
		Make Sure the door can open all the way	
		·	
ZZ-ZZ-ZZZ		/	1
10-EB-13X2			1 1
		Cabinet Light(s): "DAY BRIGHT" LED White Lights, X-Series Units	
05-IL-13S1		Cab Interior Lights Switching: Thru MultiPlex System	1
10-EB-13X3		Cabinet Light(s) : "DAY BRIGHT" LED White Light, EACH Additional,	4
		X-Series Unit	
10-EB-13YZ	S	/ Daybright Light Location Cabinet L, I-2, D2, & J	1
10-EB-13Y0	O		1 1
		Daybright location 1	
10-EB-13YB		Daybright Light Location Cabinet B	1
10-EB-13Y1		Daybright location 2	1
10-EB-13YD	S	Daybright Light Location Cabinet D1	1
10-EB-13Y2		Daybright location 3	1
10-EB-13YO		, , , ,	
		Daybright Light Location Cabinet O1	
10-EB-13Y3		Daybright location 4	1
10-EB-13YK		Daybright Light Location Cabinet K	1
10-EB-13Y4		Daybright location 5	1
10-EB-13YI	S	/ Daybright Light Location Cabinet F	1
10-EC-30A0	S		1
10-LC-30A0	3		,
		Approx. Dimensions: 17.5" W x 10" H x 29" D	
10-EC-3200	S	Location: Above Transverse O2 Below Talk Thru Window	1
10-EC-A100	S	Accessory Load Manager: For Cooling Cabinet or Fridge and	1
		Medixsafe	
ZZ-ZZ-ZZZ			1
		•	1 1
10-IL-02ST		Action Area Light: 12V, LED, Surface Mount	
10-IL-02L7		Light Location: Action Area	1
		UPHOLSTERY - CELL 7	1
11-00-UC01		Upholstery Color: Dark Gray (Gunmetal)	1
11-00-UC11		Center Trough Upholstery Color: Color Key to Rest of Truck	1
ZZ-ZZ-ZZZZ			1
			1
11-00-UT0X		Uph Joint Type: "PREMIUM " Upgrade	1
11-SB-0005		Squad Bench seat cushion cut-outs: None	1
11-SB-0003		Post And Wheel Cups: NONE	1
11-SE-CSPD		Head Protection: Pad over CS Entry Door	1
11-SE-REPD		Head Protection: 2" Pad over Rear Access Doors, Full Width	1
I I OL KLI D		PAINT - STRIPES - DECALS	
40 DT LINOT			1 1
12-PT-UNCT		Undercoating : Per QVM Guidelines, STD	1
12-00-0001	S	Wheel Finish: Painted BLACK (6) Wheels ILOS	1
12-00-0001	X		1
12-DC-0001		Reflective Tape: On painted edges of Exterior Door Frame	1
12-DC-GA10		Compartment Finish: Diamond Plate Standard	1
	C		'
12-PT-5200	3	Primary (Over All) Color: Red and Black (See Mock Up on File)	1
		Red - Paint code#FLNA3226.	
		Black - Ford OEM black (shadow black)	
		Includes Cab over all	

07/10/2019 13	:48:	31	Page 34
PART NO	S	DESCRIPTION	QTY
		Paint Entire Type 1 Red, Paint box roof White	
		Paint cab roof black.	
		Black belt above drip - rail to extrustion	
		6" Black VINYL belt with kickup (Last Unit had Black Paint, This is a change)	
12-PT-0112	S	/ Tape: Inkjet Wrap Around SOL Install Inkjet Wrap Around SOL	1
11-YZ-0962 12-PT-5500	S	Non-Compliant toKKK-F 3.16.2 ;Exterior Safety stripe to be Orange Box Roof Oxford White (YZ) ILO Over All, Not to Wrap Around Sides Roof to be painted White on Flat Section Only Not To Wrap Around Sides.	1 1
12-PT-RF01		/ Roof Paint: Color and finish quality to be GLOSSY The roof shall be painted, with a glossy finish, but not sanded and buffed.	1
12-PT-RF01	X	Paint Belt: 6" Black With Kickup	1
12-PT-RF02	X	•	1
12-PT-1000		Pin Stripe: 1/2" Gold Leaf Double & 1/2" Gold Leaf	1
12-ZZ-ZZZZ		·	1
12-PT-1314	S	/ Outer Section: including doors NOTE: To cover rear windows in semi-transparent film matching chevron pattern.	1
12-PT-1328	S	Full Rear Chevron, includes doors Diamond Grade Reflective: Double Color Include Rear Kick Plate With Chevrons	1
12-PT-1380	S	Main Color - RED	1
12-PT-1382		Second Color - YELLOW	1
12-PT-135B	S	Inside of ALL Compt/Entry Doors:	1
		State Colors Required	
		1st color RED	
		2nd color YELLOW	
		NOTE: All compt/entry doors to have chevron. Full on compt doors and lower panel on entry doors.	
12-PT-2000	S	Lettering: Custom (See Mock UP on File) -Driver Side- 5" "BROKEN ARROW" Gold Engine Turn With Black outline, Inkjet	1
		Reflective Vinyl	
		5" "FIRE DEPARTMENT" Gold Engine Turn With Black outline, Inkjet	
		Reflective Vinyl 14x20" Patch, Inkjet Reflective Vinyl	
		3" "SQUAD 7" Black Reflective	
		Install Inkjet Wrap Around SOL with Gold border and Gold Snake and Staff.	

07/10/2019 13:	:48: S		Page 35
PART NO	3	DESCRIPTION -Passenger Side-	QTY
		5" "BROKEN ARROW" Gold Engine Turn With Black outline, Inkjet	
		Reflective Vinyl	
		5" "FIRE DEPARTMENT" Gold Engine Turn With Black outline, Inkjet	
		Reflective Vinyl 14x20" Patch, Inkjet Reflective Vinyl	
		3" "SQUAD 7" Black Reflective	
		Install Inkjet Wrap Around SOL with Gold border and Gold Snake and Staff.	
		-Front-	
		4" CUSTOM REVERSE AMBULANCE WITH SOL'S	
		5" "BAFD" Gold Engine Turn With Black outline, Inkjet Reflective Vinyl -Rear-	
		2.25" "BROKEN ARROW" Gold Engine Turn With Black outline, Inkjet	
		Reflective Vinyl	
		2.25" "FIRE DEPARTMENT" Gold Engine Turn With Black outline, Inkjet	
		Reflective Vinyl 2.25" "SQUAD 7" Gold Engine Turn With Black outline, Inkjet Reflective	
		Vinyl	
		Red Yellow Chevron window perf to match diamond grade.	
40 77 7777	C	1 Overall Height Sticker (Overall Height Dive 211)	
12-ZZ-ZZZZ 12-PT-DRIP	3	/ Overall Height Sticker (Overall Height Plus 2") Drip Rails: Bright Aluminum, De-burred and rounded corners	1
		Include a drip rail over the C/S Entry Door.	•
12-ZZ-ZZZZ			1
12-PT-DRIZ		DETAIL - SHIP LOOSE ITEMS Drip Rails: Detail-Refer to Paint Section for Drip Rail Information	1 1
12-61-01/12		Include a drip rail over the C/S Entry Door.	"
		, = see	
13-00-0200	-	*Check front of W/O under chassis section for Add' items*	1
ZZ-ZZ-ZZZZ 13-10-XM00	S	/ OEM Rear View Mirror (Remove and Ship Loose)	1
13-10-XIVIOU 13-10-XM06		Owner's Manual USB Flash Drive ship loose Additional Owner's Manual Card USB Flash Drive Ship loose EACH	1
13-10-XMAN		Conversion Owners Manual, Printed w/ Wiring Schematics, EACH	1
		Ship Loose	
		A printed version of the Owner's manual will be provided in an Company ring	
		binder.	
34-OS-1400		Amico O2 Minder, Transducer: NONE	1
13-00-0001	_		1
13-00-0004	S	Decal: "CLEARANCE HEIGHT" On Windshield (Overall Height Plus 2")	1
		(Overall height Flus 2)	
<i>77-77-777</i>		/	1
13-GL-0000		Glove Holder, Ship Loose: NONE	1
ZZ-ZZ-ZZZZ			1
13-KK-KDEC 13-RF-STAR		Decals: KKK / DOT Pkg, Blue/White reflective - Ship Loose / Decals: Install 32" Star of Life on Mod Roof	1
ZZ-ZZ-ZZZZ		Decais. Instail 32 Star of Life of Wood Roof	1
13-KK-NSSB		Decals: NO SMOKING & SEATBELT, installed, cab & pt. area.	1
ZZ-ZZ-ZZZZ			1

07/10/2019 13:48:		Page 36
PART NO S	DESCRIPTION	QTY
13-KL-0010	AEV Logos: Installed on unit per AEV standard locations	1
13-MH-05ST	Fire Extinguisher, 5 pound, shipped loose, Std	1
10 1111 1001	Part No: EXTI-01	•
	Taltino. EXTI-OT	
40 141 4040		
13-MH-10A0	Reflector Pkg: Body - 2ea, Side Fr Amber, Side Rear Red, Rear Red	1
ZZ-ZZ-ZZZZ		1
13-SB-0003	Post And Wheel Cups: NONE	1
13-SP-TIRE	Spare Tire: Refer To Chassis Section to See if it gets one.	1
34-XX-1521	Primary Cot: Stryker, 6506 Power Pro XT	1
34-XX-1321		,
	This item priced for Domestic USA usage only.	
	Standard Features included :	
	1. English Language (#6506-600-000)	
	2. In Service DVD (#6506-001-102)	
	3. (2) year complete (1-yr soft goods) warranty	
	4. (3) year X-frame warranty	
	5. (3) year limited Powertrain warranty	
34-XX-1533	Smart Chip Charger System (STD)	1
	Includes (1) SMRT Charger, (2) SMRT PAKs and (1) Selected Power Cord	_
	morados (1) civil citargor, (2) civil citara (1) colocida i civor cora	
24 VV 4526	1 42V DC Cable (Automotive) #6500 027 000	4
34-XX-1536	12V DC Cable (Automotive) #6500-027-000	1
	Includes (1) SMRT Charger, (2) SMRT PAKs and (1) Selected Power Cord	
34-XX-153X	Stryker XPS side rail 6506-040-000 (note part number change)	1
	XPS is an expandable patient surface that is mounted to the sides of the cot.	
	requires XPS bolster mattress	
	requires At 6 boister mattress	
04.1/1/ 0005		
34-XX-2005	Cot Option: Stryker, IV Poles 3 Stage, Pt Right & Left Side	1
34-XX-2019	Cot STD : Stryker, Trendelenburg Feature	1
34-XX-2131	Cot Option: Stryker, Equipment Hook 6500-147-000	1
	Purchasing to include on order of patient cot configuration. 6500-147-000	
24 VV 2454	Cot Ontion: Struker Booksted Booksted Bouch 6500 120 000	4
34-XX-2151	Cot Option: Stryker, Pocketed Backrest Pouch 6500-130-000	1
34-XX-2152	Cot Option: Stryker, Head End Storage Flat 6500-128-000	1
34-XX-2250	Cot Option: Stryker, Mattress, Bolster	1
34-XX-2403	Cot Option: Stryker, Head End O2 Bottle Holder 6085-0460-00	1
	purchasing to include on patient cot at time of order.	
	1	
34-XX-2450	Cot Option: Stryker, Restraint Package, Standard	1
34-XX-2454	Cot Option:Stryker X- Restraint Package J3027 Tested SAE	1
	Compliant 6500-001-43	
	Stryker certified X restraints with the tested cot meet dynamic crash standards	
	AS/NZS-4535, BS EN-1789 and SAE J3027 as presented by Stryker EMS with a	
	Compliant Cot mount system.	
	Compliant Got mount dystom.	
24 VV 2500	L Cat Ontion Stoor Look SEGS 029 000	
34-XX-2590	Cot Option Steer-Lock 6506-038-000	1
	Purchasing activity, Required to be ordered with Cot.	
	Stryker part 6506-038-000	
	Only available for power cot with this option code	
34-XX-2750	Cot Option: Stryker Hook, "J" Hook	1
		1
34-XX-3060	Cot Option: Stryker PowerLoad Cot Compatibility 6506-127-000	T

07/10/2019 13:4		Page 37
	S DESCRIPTION	QTY
34-XX-03A5	Stryker: Stair-PRO Chair with tread system 6252	1
	Includes: This item priced for Domestic USA usage only.	
	This item priced for bornestic bort asage only.	
34-XX-03AA	Stair Chair Option : (1) Piece Vinyl Seat #6250-020-000	1
34-XX-03AC	Stair Chair Option : Locking Rear Lift Handles #6250-031-000	1
34-XX-03AE	Stair Chair Option : Extended Length Lower Lift Handles	1
34-XX-03AG	#6250-025-000 Stair Chair Option : Non Absorbent Restraint Set #6250-161-000	1
34-XX-03AG	(Metal Lap Buckles)	'
34-XX-03AJ	Stair Chair Option : In-Service Video (DVD) #6250-001-162	1
	(Metal Lap Buckles)	
34-XX-03B0 34-XX-5000	6252 Stair Chair Option : No IV Pole Clip #6252-024-000	1
	Indemnification Statement	1
	The purchaser agrees to defend, indemnify and hold REV harmless from any claims, costs (including actual attorneys' fees), damages and liabilities caused in	
	whole or in part by any alteration or modification of, or changes or additions to the	
	Purchased Products OR use of product for purposes it was not designed or	
	intended for.	
	GENERAL CONVERSION	1
36-AA-0001	WORK ORDER STATUS:	1
00-90-0304	For Proposal	1
36-CO-0001 36-CO-0010	CHANGE ORDER POLICY After Confirmation, Prior to Engineering	1
30-00-0010	Change Fee is \$50 Per Item, Plus Option Cost	,
36-CO-0020	After Confirmation, After Engineering Start	1
	Change Fee is \$100 Per Item, Plus 150% Option Cost	
	Special order or non-stocking parts that are purchased for any ordered option that is subsequently deleted will be charged for and shipped loose.	
	is subsequently deleted will be charged for and shipped loose.	
36-CO-0030	After Confirmation, After Production Start	1
30-00-0030	Change Fee is \$100 Per Item, Plus 200% Option Cost	,
	Special order or non-stocking parts that are purchased for any ordered option that	
	is subsequently deleted will be charged for and shipped loose.	
<u> </u>	!	1
36-EO-0001	<i> </i> ***END OF ORDER***	1
		1